



APPLICATION INSTRUCTIONS

JOBST® FarrowWrap® - Lower extremities



BEFORE YOU START:

- Before using compression therapy, always contact a healthcare professional. This application description does not replace instruction by an experienced expert.
- When using a JOBST FarrowWrap foot and leg section, it is recommended that you wear a knee-length JOBST Farrow Liner (understocking) or when using a JOBST FarrowWrap foot, leg, thigh and knee section wearing of a thigh-high JOBST Farrow liner is recommended.
- This application descriptions can be used for JOBST FarrowWrap LITE and JOBST FarrowWrap CLASSIC. In the pictures, JOBST FarrowWrap LITE is shown.
- For JOBST FarrowWrap LITE, BASIC and STRONG place the bands with a moderate stretch and for JOBST FarrowWrap CLASSIC the bands with a slight stretch. Please note you will see that the JOBST FarrowWrap CLASSIC thigh part includes JOBST FarrowWrap STRONG knee section with a moderate stretch.
- All jewelry should be removed from the application area.
- Open wounds in the area of application must be covered with a suitable wound dressing before application.

ORDER: First place the matching JOBST Farrow Liner on. Then start with the JOBST FarrowWrap foot part and next the JOBST FarrowWrap leg part. Before the JOBST FarrowWrap knee section, apply the thigh section.

JOBST Farrow Liner Application



Gently pull the liner over your toes and up your leg. Smooth any wrinkles. (With knee-length liners (1a), give your toes room for movement by pulling the fabric slightly away from the toes).


JOBST FarrowWrap Footpiece Application



Place your heel into the heel pocket. Gently stretch with medium tension and affix the foot flaps.



Gently stretch the ankle bands over the top of the foot to the opposite side and affix, ankle should be at 90 degrees flexion.



JOBST FarrowWrap Legpiece Application



Lay the wrap behind your leg with the label side facing up.

Tip: As demonstrated in the image, for easier handling you can start by closing the top band below your knee first.



Grab each side of the band closest to your ankle and pull forward wrapping the band around the lowest part of the leg by applying a gentle stretch and affix with the VELCRO® band. If wearing a footpiece, this band should overlap the footpiece. Apply the remaining bands in the same manner.



Repeat with each band moving up the leg. When properly donned, the white label will be located just below the knee crease.

JOBST FarrowWrap Thighpiece Application



Take the thigh piece and place it underneath the thigh so that it rests just above the knee with the label centered on the back of the upper thigh. Wrap the lower thighpiece band around the lower thigh, just above the knee and by applying a gentle stretch affix it with the VELCRO®.



Apply the remaining band(s) in the same manner.

JOBST FarrowWrap Kneepiece Application



Place the middle band of the kneepiece underneath the knee, with or without the foam pad. Wrap the middle band around the knee by applying a moderate stretch.



Place the remaining bands in the same manner. The middle band should be visible through a gap in the middle of your kneecap.

NOTES:

- JOBST FarrowWrap must be regularly readjusted if necessary, to ensure the correct fit.
- To ensure that the Velcro fasteners do not get caught when putting them on, they can be rolled onto the back, rolled up and fixed in place.
- If desired, fold the protruding liner over the upper edge of the JOBST FarrowWrap leg or thigh section.



TIPS FOR CORRECT APPLICATION



Washing instructions

Close the Velcro fasteners before washing to prevent them from snagging or getting caught. Always use a laundry net to protect the product.



Covering wounds

Before applying JOBST FarrowWrap, you should ensure that all wounds are thoroughly cleaned. All open wounds present in the area to be treated must be covered with a suitable wound dressing before applying the compression garment.



Wear suitable footwear

Safe footwear is important to reduce the risk of falls and to protect the JOBST FarrowWrap medical adaptive compression systems from premature wear and tear.



Storage

To ensure a longer service life of the compression garment, JOBST FarrowWrap should be stored properly. We recommend folding the Velcro fasteners over to the back and closing them, to avoid snagging.

For more information www.jobstcanada.com