# TENA SmartCare Change Indicator™/MC

|    |   |        |
|----|---|--------|
| en | Family caregivers Instructions for use                | 2–35   |
| es | Familiar cuidador Instrucciones de uso                | 36–69  |
| fr | Instructions d'utilisation pour les aidants familiaux | 70–104 |

This guide is intended for family caregivers using TENA SmartCare Change Indicator<sup>TM/MC</sup> in a home environment.

TENA SmartCare Change Indicator is a digital health technology solution that helps family caregivers decide when to change absorbent products without the need for frequent manual checks.

The Change Indicator is attached to the outside of the absorbent product and tracks the urine saturation levels within. Family caregivers then receive information about the different saturation levels on their smartphones via the TENA SmartCare Family Care app. This information is then used to decide whether the absorbent product needs to be changed or not. The TENA SmartCare Change Indicator system tracks urine saturation, but does not detect fecal episodes.

**TENA SmartCare Change Indicator is a system consisting of:**

- An app named TENA SmartCare Family Care that can be installed on one or more smartphones used in a family.
- A reusable Change Indicator that consists of two parts: a sensor strip and a transmitter.
- A gateway that relays data from the Change Indicator to a secured server.

**In this guide, you will find instructions for how to use:**

The TENA SmartCare Change Indicator system including TENA SmartCare transmitter, TENA SmartCare gateway, TENA SmartCare sensor strip and TENA SmartCare Family Care app.

The TENA SmartCare Change Indicator system is designed to be intuitive, with family caregivers in mind. It does not require special training.

**Symbols used in this guide**

**WARNING** Care shall be taken to ensure safety and efficiency. Failure to do so could cause damage to equipment or lead to the receipt of incorrect data.

**NOTE** Important information or advice/ tips.

**CHECK** Important information for the system to function.

|  |              |
|--|--------------|
| <b>Important information</b>   | <b>4</b>     |
| <b>Warnings</b>  | <b>5</b>     |
| <b>Included in the TENA Smartcare Change Indicator™/MC starter kit</b> | <b>7</b>     |
| <b>Download the app and set up the account</b>                         | <b>8</b>     |
| <b>Explanation of the app during use</b>                               | <b>9</b>     |
| <b>Set up the gateway</b>  | <b>10</b>    |
| <b>Register the transmitter code</b>                                   | <b>11</b>    |
| <b>Install the battery &amp; assemble the Change Indicator</b>         | <b>12</b>    |
| <b>Attach the Change Indicator</b>                                     | <b>13</b>    |
| <b>Use the Change Indicator</b>  | <b>14</b>    |
| <b>Remove the Change Indicator</b>                                     | <b>15</b>    |
| <b>Intended use</b>  | <b>16</b>    |
| <b>List of supported TENA products</b>                                 | <b>17</b>    |
| <b>Attachment instructions for supported TENA products</b>             | <b>18</b>    |
| TENA Underwear   | 18           |
| TENA Briefs  | 19           |
| TENA Belted briefs   | 20           |
| <b>Replace the sensor strip &amp; the battery</b>                      | <b>21</b>    |
| <b>Cleaning</b>  | <b>22</b>    |
| <b>Gateway lights &amp; buttons</b>                                    | <b>23</b>    |
| <b>Transmitter lights &amp; buttons</b>                                | <b>24</b>    |
| <b>Troubleshooting</b>   | <b>25–29</b> |
| <b>Technical descriptions</b>  | <b>30</b>    |
| Operating characteristics  | 30–31        |
| Electromagnetic compability (EMC)                                      | 32           |
| Lifetime   | 32           |
| Disposable parts   | 33           |
| Labelling symbols  | 34           |
| <b>Assortment overview</b>   | <b>105</b>   |

- Ensure you have read and understood these instructions before using the product. Save these instructions for future reference.
- If you encounter any problems, equipment failure or if any other incidents occur, report them to the manufacturer or TENA SmartCare support.
- The Change Indicator is intended for use with TENA incontinence products.
- The adhesive strip for attaching the gateway will not adhere well to dirty, dusty or oily surfaces and may not adhere well to vinyl surfaces. Do not use tape on wallpaper, bare wood or on porous or painted surfaces.
- The Change Indicator is not intended to be used in direct contact with the skin.
- Do not use gateways with unapproved power supplies.

- Keep the Change Indicator and batteries away from children and pets.
- Only a CR2032 battery should be used in the Change Indicator. Using the wrong type of battery, including rechargeable batteries, can damage the Change Indicator.
- The Change Indicator should be attached to the outside of the incontinence product.
- Clean the Change Indicator when fully assembled to prevent water from entering the device. Do not submerge in water or clean in a washing machine.
- Do not swing the Change Indicator in the direction of people, pets or nearby objects. Be aware of the risk of strangulation with the sensor strip.
- Seek immediate medical care if the battery is swallowed.
- Do not attempt to repair or modify the Change Indicator or a gateway by yourself.
- If any part of the Change Indicator system appears to be damaged or broken, it should not be used.
- Do not expose the battery to high temperatures, disassemble it, damage it mechanically or put battery into the fire as this can cause an explosion or a fire.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:


- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.


This device contains license-exempt transmitter(s)/ receiver(s) that comply with Innovation, Science and Economic Development Canada's license-exempt RSS(s). Operation is subject to the following two conditions:

- This device may not cause interference.
- This device must accept any interference, including interference that may cause undesired operation of the device.

This equipment complies with the ICES RF radiation exposure limits set forth for an uncontrolled environment.

**You also need:**

**A SUPPORTED TENA INCONTINENCE PRODUCT**  
(see supported TENA products on p. 17 of this guide or in the TENA SmartCare Family Care app).


**A SMARTPHONE**  
with TENA SmartCare Family Care app, iOS 11/ Android 7 or later.

**CONSUMABLES:**

Sensor strip 1 (short), sensor strip 2 (medium), sensor strip 3 (long), battery, TENA incontinence product.


**1** **DOWNLOAD THE TENA SMARTCARE FAMILY CARE APP**

Download the TENA SmartCare Family Care app from the App Store or Google Play and then register your account.


**2** **INVITE OTHER CAREGIVERS**

You can create a care team by inviting others who care for your loved one. This will give them access to the TENA SmartCare Family Care app.

**3** **REGISTER THE TRANSMITTER**

Follow the setup instructions in the TENA SmartCare Family Care app. The transmitter code is located on the back of the transmitter and on the packaging.

**1 DRY**  
- no action needed

The absorbent core of the TENA product is dry, and the surface feels dry to the skin.

**2 SOME URINE**  
- no action needed

There is some urine in the absorbent core of the TENA product. However, the surface still feels dry to the skin. If the absorbent product has been on for an extended period and the user has sensitive skin, you may want to consider changing it.

**3 WET**  
- consider change

There is a large amount of urine in the absorbent core of the TENA product. The surface probably still feels dry to the skin, but you should consider a change - especially if the absorbent product has been on for an extended period and/ or the user has sensitive skin.



**NOTE**

The gateway should be placed within 33 ft/ 10 m of the Change Indicator and connected to an electrical power outlet.

**Placement**


**1** Find the transmitter code, which is located on the back of the transmitter and on the packaging.

**2** Enter the transmitter code in the TENA SmartCare Family Care app to register.

# en INSTALL THE BATTERY & ASSEMBLE THE CHANGE INDICATOR

## Install the battery

Light blinks  
5 times = OK



See page 24  
Transmitter lights  
& buttons

**WARNING** Only CR2032 batteries shall be used for TENA SmartCare Change Indicator.

## Assemble the Change Indicator (transmitter + sensor strip)


Steady light for a  
few seconds = OK


See page 24  
Transmitter lights  
& buttons

See supported TENA absorbent products on p. 17 of this guide or in the TENA SmartCare Family Care app.


**1 Find the front**


**2 Keep folded**

**3 Align with the top of the thicker absorbent area**


**4 Place tightly**


**NOTE**

Always attach the Change Indicator before you unfold the incontinence product. Start by placing the Change Indicator at the front of the incontinence product, on the thicker absorbent area, ensuring it is straight.

Align the transmitter with the top of the thicker absorbent area of the incontinence product. Press firmly to ensure the Change Indicator has full contact with the incontinence product. If the user is very mobile, consider applying a fixation pant or ordinary underwear on top of the Change Indicator. This will keep it comfortably in place.

1

Press both buttons simultaneously for two seconds at every change. This records the change of incontinence product in the TENA SmartCare Family Care app.

Steady light for a few seconds:  
Continue



Blinking light:  
Change battery, see page 21

2

Apply the TENA absorbent product with the Change Indicator on your loved one.

3

For explanations of app notifications, see page 9.

**DRY**

- no action needed

**SOME URINE**


- no action needed


**WET**

- consider change

1

Start at the top. Grip the tab, then gently peel straight downwards.


Don't grab and pull the transmitter.

**NOTE**

The Change Indicator should be removed while the user is still wearing the incontinence product. Gently pull the Change Indicator down from the top, then roll the user onto his/ her side before completing the removal.

2

Always make sure the Change Indicator has been removed before disposing of the TENA incontinence product.

**Intended use:**

The TENA SmartCare Change Indicator™/MC is an accessory to absorbent incontinence products, intended for use on individual(s) suffering from urinary incontinence in a home or professional environment who are dependent on one or more caregivers to change the absorbent incontinence products. The TENA SmartCare Change Indicator estimates the degree of urine saturation in the absorbent incontinence product and notifies the caregiver(s). This facilitates the caregiver decision regarding when to change the absorbent incontinence product.

**Standard operation:**

- If the Change Indicator is broken, it should be replaced, not repaired.
- The Change Indicator should not be in contact with the user's skin during normal use. It should always be attached to the outside of the incontinence product.
- The Change Indicator is operated by the caregiver of the incontinent user. The user wears the Change Indicator in combination with a TENA incontinence product.
- Attaching the Change Indicator to the incontinence product is considered a primary operating function.
- Only the caregiver is intended to interact with the Change Indicator.
- To turn the transmitter off, remove the battery.
- The Change Indicator is intended for persons with healthy skin.

**Advanced operation:**

- Changing the battery and sensor strip are operations that should be carried out infrequently.

Using the TENA SmartCare Change Indicator according to its intended use does not add any risk compared to not using it.

The medical device complies with emission and immunity standards and regulations. It has been tested and fulfills requirements for home healthcare environments.

**CHECK** Ensure you have the correct sensor strip length for your TENA product. Only products with a textile-like outer lining can be used with the Change Indicator.

| TENA®<br>Absorbent<br>Products                                 | TENA®<br>Productos<br>absorbentes                                      | Produits<br>absorbants<br>TENA®  | S   | M   | R | L    | XL        | 2XL         | 3XL |
|--|--|--|-----|-----|---|------|-----------|-------------|-----|
|  |  |  | S/M | M/R |   | L/XL |           |             |     |
|  |  |  | P   | M   |   | G    | XG,<br>TG | 2XG,<br>2TG |     |
| <b>Briefs</b>  | <b>Pañales</b>   | <b>Culottes</b>  |     |     |   |      |           |             |     |
| TENA Stretch™/MC: Ultra / Super                                |  |  |     | 2   |   | 3    |           | 3           | 3   |
| TENA®: Ultra / Super   |  |  |     | 2   | 2 | 3    | 3         |             |     |
| TENA Complete™/MC /<br>TENA Complete + Care™/MC                |  |  |     | 2   |   | 2    | 2         |             |     |
| TENA Complete Ultra™/MC /<br>TENA Complete + Care Ultra™/MC    |  |  |     | 2   |   | 2    | 2         |             |     |
| TENA® Small<br>Brief   | TENA® Pañal<br>Pequeño   | TENA® Culottes<br>Petites  | 1   |     |   |      |           |             |     |
| TENA Bariatric Ultra™/MC                                       |  |  |     |     |   |      | 3         |             |     |
| <b>Underwear</b>   | <b>Ropa interior</b>   | <b>Sous-vêtement</b>   |     |     |   |      |           |             |     |
| TENA® Underwear:<br>Plus / Extra /<br>Overnight Super          | TENA® Ropa<br>Interior:<br>Plus / Extra /<br>Overnight™ Super          | TENA®<br>Sous-vêtement<br>Absorbant:<br>Plus / Extra /<br>Overnight™ Super | 1   |     |   |      |           |             |     |
| TENA ProSkin™/MC<br>Protective<br>Underwear for<br>men / women | TENA ProSkin™<br>Ropa Interior<br>Protectora para<br>hombres / mujeres | TENA ProSkin™<br>Sous-vêtement<br>pour hommes /<br>femmes                  | 1   |     |   |      |           |             |     |
| <b>TENA ProSkin™/MC Belted brief</b> : Super, Maxi             |  |  | 1   | 2   |   | 3    | 3         |             |     |


## TENA Underwear

The TENA SmartCare Change Indicator™/MC should always be attached when the incontinence product is folded and before it has been applied to the user.


### NOTE

Ensure the Change Indicator is straight. Smooth any sensor strip wrinkles, then press firmly to ensure the Change Indicator has full contact with the TENA absorbent product.

1

**Find the front**


2

**Keep folded**


3

**Align with the top of the thicker absorbent area**

4

**Place tightly**


5

Every time you change an incontinence product, press both buttons simultaneously for two seconds to reset the transmitter and register the change in the app.

6

Apply the absorbent product with the Change Indicator to your loved one.


## TENA Briefs

The TENA SmartCare Change Indicator™/MC should always be attached when the incontinence product is folded and before it has been applied to the user.

### NOTE


Ensure the Change Indicator is straight. Smooth any sensor strip wrinkles, then press firmly to ensure the Change Indicator has full contact with the TENA absorbent product.

1

### Find the front

FRONT


2

### Keep folded


3

### Align with the top of the thicker absorbent area


4

### Place tightly


5

Every time you change an incontinence product, press both buttons simultaneously for two seconds to reset the transmitter and register the change in the app.

6

Apply the absorbent product with the Change Indicator to your loved one.


# TENA Belted Briefs

The TENA SmartCare Change Indicator™/MC should always be attached when the incontinence product is folded and before it has been applied to the user.

### NOTE

Ensure the Change Indicator is straight. Smooth any sensor strip wrinkles, then press firmly to ensure the Change Indicator has full contact with the TENA absorbent product.

**1 Find the front**

Find the front based on the application of the product.

**2 Keep folded**

**3 Align with the top of the thicker absorbent area**


**4 Place tightly**

Every time you change an incontinence product, press both buttons simultaneously for two seconds to reset the transmitter and register the change in the app.

Apply the absorbent product with the Change Indicator to your loved one.


### Sensor strip replacement

The sensor strip is designed to withstand approximately 12 weeks of continuous use. Sometimes, however, the hooks on the back of the sensor strip can become covered with fibers or lint. This prevents the hooks from attaching to the incontinence product. If the fibers or lint cannot be removed, the sensor strip should be replaced. The sensor strip should also be replaced if it becomes soiled with feces and cannot be cleaned sufficiently. For cleaning advice, see page 22 of this guide.

### Battery replacement

If the LED light starts to blink, either the battery needs to be replaced or the transmitter is not properly attached to the sensor strip. The CR2032 battery will last for up to 12 weeks. However, the battery can be replaced whenever needed.


#### CHECK

If you have difficulty removing the battery, first press the locking tab, then tap the transmitter in your hand until the battery falls out. The transmitter is designed to withstand attachment and removal from a sensor strip four times per year over two years. With more frequent attachments/ removals, the system may not work as intended.

### Change Indicator

The transmitter and sensor strip should never be cleaned separately – only clean the Change Indicator when it is fully assembled. Clean the Change Indicator using a cloth or a disinfectant when required.


We recommend using:

- Pre-soaked, alcohol-based wipes without additives.
- Alcohol-based disinfection foam for the hooks on the back side of the sensor strip.


The assembled Change Indicator can be rinsed under a tap.

Any fibers or lint attached to the hooks on the back of the sensor strip should be removed to ensure a secure fit between the Change Indicator and the incontinence product. Fibers can be removed using the free end of the sensor strip.

### Gateway

Clean the gateway with a dry cloth when needed.


#### **WARNING**

Clean the Change Indicator when fully assembled to prevent water from entering the device.


Do not clean in the washing machine.

Do not submerge in water.

### Gateway installation

After attaching the gateway to a power outlet:


**NOTE**

The LED light can be difficult to see if exposed to a strong sunlight, for example. Try to block the sunlight to improve readability.

| COLOUR LED LIGHT  | WHAT IT MEANS   |
|---|---|
|   | Blinking blue light:<br>Wait. The system is connecting.<br>Allow up to 60 seconds for the light to start blinking.  |
|  | Steady green light:<br>Everything is OK. Allow up to 60 seconds for the light to turn green.  |
|  | No light:<br>Green light disappears after a while and the system goes into dark mode. Push the button to see if it has power and network connection. If no light appears, check the power connection. |
|  | Steady red light:<br>Error. See p. 25–29 for more information.  |
|  | Blinking blue-green-red light:<br>pull out the power. Wait 30 sec. Insert the power. Wait 2 min.<br>If the error persists, contact TENA SmartCare support.  |

**NOTE**

The LED light can be difficult to see if exposed to a strong sunlight, for example. Try to block the sunlight to improve readability.

**Battery installation**

Once the battery has been installed in the transmitter:

| LED LIGHT   | WHAT IT MEANS  |
|---|--|
|  | Blinks five times:<br>The battery is in the correct position. Attach the transmitter to the sensor strip.                          |
|  | No light: OK<br>The light disappears after a few seconds. Check for power by pressing both buttons simultaneously for two seconds. |

**Change Indicator assembly**


Once the transmitter has been attached to the sensor strip:

| LED LIGHT   | WHAT IT MEANS  |
|---|--|
|  | Steady light for a few seconds:<br>Everything is OK.   |
|  | Blinks five times:<br>The transmitter is not attached properly or the battery is low. Try re-attaching the transmitter. If it does not help, change the battery. |
|  | No light: OK<br>The light disappears after a few seconds. Check for power by pressing both buttons simultaneously for two seconds.                               |

**Applying the Change Indicator to the incontinence product**

Press both buttons simultaneously for two seconds at every change of incontinence product to register the change in the app.

| LED LIGHT   | WHAT IT MEANS  |
|---|--|
|  | Steady light for a few seconds:<br>Everything is OK.   |
|  | Blinks five times:<br>The transmitter is not attached properly or the battery is low. Try re-attaching the transmitter. If it does not help, change the battery. |
|  | No light:<br>Error. See pp. 25–29 for more information.  |

| PROBLEM                               | SOURCE  | SOLUTION   |
|---------------------------------------|---------|--|
| <b>Connection error (signal lost)</b> | Gateway | <ol style="list-style-type: none"><li>1. Check power supply. Are the power cable and power adapter both properly attached?</li><li>2. Press the gateway button. If there is a steady green light, everything is OK.</li><li>3. If there is a blinking blue light, the gateway is trying to connect.</li><li>4. If the gateway fails to connect and the light doesn't turn green after one minute, try moving the gateway to a position in the room where the network connection is stronger and wait for the green light.</li><li>5. Occasionally, the cellular network may go down. If this happens, revert to a manual changing routine until the system is up and running again.</li><li>6. If this does not help, please contact TENA SmartCare support.</li></ol> |

| PROBLEM                               | SOURCE           | SOLUTION  |
|---------------------------------------|------------------|---|
| <b>Connection error (signal lost)</b> | Change Indicator | <ol style="list-style-type: none"> <li>1. Make sure the transmitter is properly attached to the sensor strip and that the battery is properly inserted (see p. 12 of this guide).</li> <li>2. Press both buttons on the assembled Change Indicator simultaneously for two seconds.               <ol style="list-style-type: none"> <li>a. If there is no light, change the battery and re-attach the transmitter to the sensor strip. Press both buttons simultaneously. If there is a steady light for a few seconds, everything is OK.</li> <li>b. If there is a blinking light, check that the transmitter is properly attached to the sensor strip. Adjust if needed, then press both buttons simultaneously. If the light continues blinking, change the battery.</li> <li>c. If there is a steady light for a few seconds, everything is OK.</li> </ol> </li> <li>3. If this does not help, replace the sensor strip. Press both buttons simultaneously. If there is a steady light for a few seconds, everything is OK.</li> <li>4. If this does not help (no light or blinking light after correctly attaching to the sensor strip), please contact TENA SmartCare support.</li> </ol> |

| PROBLEM                               | SOURCE           | SOLUTION   |
|---------------------------------------|------------------|--|
| <b>Connection error (signal lost)</b> | System range     | <ol style="list-style-type: none"> <li>1. The Change Indicator should be in range of the gateway – usually within 33 feet/ 10 meters.</li> <li>2. Try moving the gateway to a place where it is within range of the Change Indicator.</li> <li>3. If this does not work, check the “connection error” section above.</li> </ol>  |
| <b>Change Indicator falls off</b>     | Change Indicator | <ol style="list-style-type: none"> <li>1. Make sure you are using one of the supported TENA incontinence products with a textile-like outer lining. The Change Indicator does not attach to absorbent products with a plastic outer lining.</li> <li>2. Make sure you attach the Change Indicator to a folded product, straight from its plastic packaging. Do not attach it when the user is already wearing the absorbent product.</li> <li>3. Check the hook-covered side of the sensor strip to see if it is covered with lint. Remove any lint.</li> <li>4. If this does not help, replace the sensor strip.</li> </ol> |

| PROBLEM  | SOURCE                   | SOLUTION   |
|--|--------------------------|--|
| <p><b>The system notifies me too early/ too late</b></p> | <p>Status indication</p> | <p>a) You receive the “consider change” notification when an absorbent product is still dry or has not been worn for long.</p> <p>b) You receive the “consider change” notification when an absorbent product is already too full of urine to be comfortable.</p><br><ol style="list-style-type: none"> <li>1. Make sure you are using one of the supported TENA products (see p. 17 of this guide).</li> <li>2. Make sure you are using the correct sensor strip length for your TENA product (see p. 17 of this guide).</li> <li>3. Make sure the Change Indicator is attached correctly (see pp. 18–20 of this guide).</li> <li>4. Check the settings in the TENA SmartCare Family Care app. Is the proper absorbent product selected?</li> <li>5. If this does not help, please contact TENA SmartCare support.</li> </ol> |

| PROBLEM   | SOURCE            | SOLUTION  |
|---|-------------------|---|
| <b>I do not receive any notifications</b>   | Status indication | <ol style="list-style-type: none"><li>1. Check “connection error” in the troubleshooting section of this guide.</li><li>2. Check that you have entered the right transmitter code in the TENA SmartCare Family Care app.</li><li>3. If this does not help, please contact TENA SmartCare support.</li></ol> |
| <b>The incontinence product was filled with feces, but the system did not notify me</b> | Feces             | The system does not detect feces. In case of frequent fecal episodes, please revert to a manual changing routine.   |

## Operating characteristics

The Change Indicator is classified as an internally powered device for continuous operation.

## Processing units

Transmitter & gateway

| CHARACTERISTICS            | VALUE/ RANGE  |
|----------------------------|---|
| Wireless technology        | RF communication, BLE   |
| RF frequency               | 2.4 GHz   |
| Max. output level (ERP)    | 0 dBm   |
| Electrical emissions       | Class B, level 47 dB mV/m   |
| Operating voltage          | 2.7 – 3.3 V   |
| Power consumption          | Avg. 135 uW<br>Peak 75 mW   |
| Remote communication range | At least 10 meters/ 33 feet in open environment, without obstructions |

## Battery module

| TYPE   | CAPACITY/ VOLTAGE      |
|--------|------------------------|
| CR2032 | 220 mAh/ 3.0 V nominal |

This value is with the recommended battery: CR2032

## Product component dimensions (mm/ inch)

| ITEM                   | LENGTH                       | WIDTH          | HEIGHT        | DIAMETER     |
|------------------------|------------------------------|----------------|---------------|--------------|
| Transmitter            | 51.8 mm/ 2.04"               | 34.5 mm/ 1.36" | 13 mm/ 0.51"  | N/A          |
| Sensor strip (1, 2, 3) | 510 – 690 mm/ 20.08 – 27.16" | 27.5 mm/ 1.08" | 7.6 mm/ 0.30" | N/A          |
| Battery                | N/A                          | N/A            | 3.2 mm/ 0.13" | 20 mm/ 0.79" |

## Product weight (g/ ounce)

| COMPONENT              | WEIGHT                        |
|------------------------|-------------------------------|
| Transmitter            | 12.6 g/ 0.44 oz               |
| Sensor strip (1, 2, 3) | 12.6 – 15.9 g/ 0.44 – 0.56 oz |
| Battery                | 3.0 g/ 0.11 oz                |

## Environmental conditions

| CONDITION                           | MINIMUM       | MAXIMUM        |
|-------------------------------------|---------------|----------------|
| Storage & transport temp.           | +10 °C/ 50 °F | +40 °C/ 104 °F |
| Storage & transport humidity        | 65% Rh        | 90% Rh         |
| Operating temp.                     | +15 °C/ 59 °F | +40 °C/ 104 °F |
| Operating humidity (non-condensing) | 15% Rh        | 90% Rh         |
| Pressure                            | 700 hPa       | 1060 hPa       |

## Electromagnetic compatibility (EMC)

RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm/ 12" to the Change Indicator. Otherwise, degradation of the performance of this equipment could occur.

Interference may occur in the vicinity of equipment marked with this symbol.

## Lifetime

- Read best before date on accompanying packaging.
- Expected lifetime of the transmitter is two years and the sensor strip is 12 weeks, after the first use.
- Expected shelf life (from the manufacturing date) for the transmitter is three years and for the sensor strip is three years.
- The transmitter is designed to withstand attachment and removal from a sensor strip four times per year over two years. With more frequent attachments/ removals, the system may not work as intended.

## Disposable parts

- Clean the Change Indicator before disposal.
- Follow local regulations for disposal of the Change Indicator, the gateway and the batteries.
- This declaration relates exclusively to the radio equipment in the state in which it was sold and excludes components that are added and/ or operations carried out subsequently by the final user.

### WEEE mark:

- This product has been designed with recycling, reuse and waste management in consideration. Essity will provide the necessary information to recycling facilities and government agencies as required by WEEE Directive. For more information, please contact the manufacturer at [www.tena.us](http://www.tena.us).
- The product will be labeled with the WEEE mark as shown on the side, in accordance with European Standard EN 50419.
- The proper disposal of WEEE products will conserve natural resources (through reuse and recycling) and protect environmental and human health through proper disposal.
- WEEE-labeled products must be kept separate and returned to a designated collection point for dismantling, recycling and/ or reuse.
- The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local government, your commercial, institutional and/ or municipal waste disposal service or the location where you purchased the product.

## Labeling symbols

The following symbols may appear on your TENA SmartCare Change Indicator<sup>TM/MC</sup>, accessories and/ or packaging:

Manufacturing date

**FCC ID:** FCC ID number

Manufacturer

**IC:** IC ID number

Instructions for use

**UDI:** Unique Device Information

Not MR safe




WEEE mark  
Dispose of electrical components in accordance with your local regulations

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Essity Hygiene and Health AB
Essity HMS North America Inc.