TENA®

smartcare

BLADDER SENSOR

INSTRUCTIONS FOR USE

Version
IFU EN 2-0

Date
2025-01-13

Contact Information
customer Service
www.tena.com

CONTENT en

1. Product description 4
2. Intended users 9
3. Intended use and medical claims 10
4. Instructions. 11
5. Service life 35
6. Contraindications 36
7. Precautionary information. 37
8. Specifications. 41
9. Symbols. 46
10. Transport Conditions 50
11. Manufacturer 52
12. EMC Report 53
13. Applicable norms for the Bladder Sensor 63

CE
1639

The TENA SmartCare Bladder Sensor is designed to assist children and adults who struggle to recognize when their bladder is full.

It provides real-time information on the fullness of the bladder and informs the user when it's time to go to the toilet.

The device gently vibrates on the belly and sends a notification to the user's smartphone and/or Apple Watch via Bluetooth®.

To stick the device onto the belly, the user will use an adhesive and apply gel between the device and the skin.

1.1 INCLUDED IN THE PACKAGING (A-F)

A. TENA SmartCare Bladder Sensor with ultrasound sensor


Latest software version included in the app

Manufacturer
Novioscan B.V.
Transistorweg 5,
6534 AT Nijmegen
The Netherlands

Reference/Type
Product code:
08719326106298

B. USB charger + blade(s)

Only supplied charger should be used

Manufacturer
GlobTek, Inc.
www.globtek.com
186 Veterans Drive,
Northvale, NJ 07647
T: (201) 784-1000

Reference/Type
GTM46101-1005-USB

C. USB charger cable

Only supplied charger cable should be used.

Manufacturer

BENELUX rep.:
Allekabels BV, Postbus 1077
6460 BB Kerkrade
The Netherlands
T: +31 (0)45 5717001

Reference/Type

USB TYPE-C CABLE,
20CM, USB 3.1 TO
USB 3.0A

D. App[1]

To install the app, you need a Smartphone and/or Apple Watch with:

- Bluetooth® Low Energy
- iOS 16 or higher
- WatchOS 9.3 or higher
- Android 12 or higher

Manufacturer

Novioscan B.V.
Transistorweg 5
6534 AT Nijmegen
The Netherlands


Reference/Type

Download the latest version in the App Store.

[1] *Inseparable part of the medical device. TENA SmartCare Bladder Sensor.*

E. Quick guide manual

Quick guide workflow for users.

Manufacturer

Novioscan B.V.
Transistorweg 5,
6534 AT Nijmegen
The Netherlands

Reference/Type

Version 1-0

F. Instructions for use

Complete instructions for use.

Manufacturer

Novioscan B.V.
Transistorweg 5,
6534 AT Nijmegen
The Netherlands

Reference/Type

Version 2-0

1.2 ADDITIONALLY REQUIRED FOR USE (G-H)

Only Novioscan adhesive patches may be used.

Manufacturer

Novioscan B.V.
Transistorweg 5
6534 AT Nijmegen
The Netherlands

Reference/Type

Product code:
08719327757208

H. TENA SmartCare Ultrasound Gel

Only GELLO® for Ultrasonic Use may be used.

Manufacturer

GELLO® GmbH Geltechnik
Harmate 28, 48683 Ahaus
Germany
T: +49 2561/9802-0

Reference/Type

GELLO® ultrasound gel
Art.-No: 04251765105035

The TENA SmartCare Bladder Sensor is designed to assist children (aged 6 and above) and adults (with a BMI range of 18.5 kg/m^2 to 25 kg/m^2) who experience urinary incontinence during the day and/or night.

It can be used in both home and healthcare settings, such as hospitals and nursing homes.

The device informs users when it is time to use the toilet by vibrating gently on the abdomen and sending a notification to the app via Bluetooth® on a smartphone and/or smartwatch.


Caregivers can also use the app to remind or assist the user to use the toilet.

The TENA SmartCare Bladder Sensor is a Class IIa device as specified in the EU Medical Device Regulations (2017/745).

TENA SmartCare Bladder Sensor is an on the body ultrasound device for continuous bladder monitoring for support in the prevention of urinary incontinence, in healthcare and homecare environment.

The intended clinical benefits of the Bladder Sensor are to support children and adults who suffer from urinary incontinence during day and/or night. Knowledge of when the bladder has reached, or is close to reaching, its maximum capacity, will support users with urinary incontinence in understanding the true sense of urgency and allow them to go to the bathroom on time and thereby reduce the number of urinary leakages. In addition, data collected by the Bladder Sensor will provide insights into the voiding behaviour of the user. Next, wearing Bladder Sensor can contribute to a positive change in the quality of life of users such as being less affected by urinary incontinence during normal daily activities, sports and/or social life, worrying less regarding smell and feeling less embarrassed.

4.1 DEVICE LAYOUT

Battery light:

The light represents the battery power.

Bluetooth® light:

The light indicates whether the device is connected by Bluetooth® to the app.

Position search light:

The light indicates whether the bladder is detected and is only active during the position check.

Sensor:

Ultrasound gel needs to be applied to this area.

BATTERY LIGHT

Green blinking	Battery power is above 15%.
Red solid	Battery power is *below* 15%. Recharge the sensor.

BATTERY LIGHT WHILE CHARGING

Green	Battery is charged.
Orange solid	Battery is recharging. Keep charging.

INDICATIONS BY POSITION SEARCH LIGHT

Red solid	The Bladder Sensor is upside down (also indicated in the app).
Red solid	No bladder is seen.
Green blinking	Bladder is seen.
Green solid	Bladder is seen longer than 3 seconds.

BLUETOOTH® LIGHT

Blue blinking	Device is connected.
Off	Device is not connected.

4.2 GETTING STARTED

WARNING: Minors must be aided by an adult during the setup and use of the sensor.

Required materials

- Bladder Sensor
- USB charger (incl. blade) and cable
- Adhesive patches
- Ultrasound gel
- Smartphone and/or Apple Watch
- Bladder Sensor app

Charge the Bladder Sensor until the battery light is green

Note: During charging, the sensor is inactive. The device is disconnected from the electrical supply by unplugging from the charger.

The sensor starts automatically after unplugging.

Install app

Go to the app store on your Smartphone to download the Bladder Sensor app.

4.2 GETTING STARTED

Open the app & activate Bluetooth®

Note: You can open the app on your Smartphone and/or Apple Watch.

Follow the onboarding steps presented in the overview of the app:

Protection against unauthorized access

Note: You can only open the App on your smartphone or smartwatch if your mobile device is secure (i.e. passcode is set). To continue without passcode, user consent is required taking full responsibility for any consequences of a potential security breach.

Step 1. Connect your sensor

A Press the Bluetooth® icon to search for the sensor.

B Select the sensor that you wish to connect to.

Note: Not able to connect the sensor to your Smartphone and/or Apple Watch? Please plug & un-plug your sensor from the charger.

The number of your sensor can be found on the label (behind [SN]):

Step 2. Bladder size setup

The sensor has a default threshold setting size 4 (175–225 mL) > see table 1. The bladder size can be adjusted later in the settings of the app (see p .17 for an overview of all threshold settings corresponding to the bladder capacity).

Table 1: Threshold settings for bladder size

BLADDER SIZE	EXPECTED BLADDER CAPACITY – mL
1	0 – 75
2	75 – 125
3	125 – 175
4	175 – 225
5	225 – 275
6	275 – 325
7	325 – 375
8	375 – 425
9	425 – 475
10	475 – 525

4.2 GETTING STARTED

Step 3. Attach your sensor

Follow the attachment guide presented in the app:

Ultrasound gel application

Apply gel to the black rectangular area on the back of your sensor and make sure the full area is covered with gel.

Position check

The next screen should lead to the position check. If not, contact customer service.

Reminder: A full bladder is required to do the position check. To get started, drink 2 glasses of water and wait for 45-60 minutes before following the instructions on the next page.

The detection of the bladder may be impacted by anatomical limitations e.g., loose skin or a protruding belly. The position check helps you to find the correct position of the sensor on the lower abdomen. The app will guide you through the position check.

Prepare position check

Hold the sensor with gel against your lower abdomen and approximately 1 cm above your pubic bone to find the right spot where the sensor can locate your bladder.

You can feel the pubic bone if you press your fingers against your lower abdomen.

Locate your bladder

Activate the position check in the app. The sensor search light will blink green when the bladder has sucessfully been detected.

Make sure your bladder is not empty.

The sensor should be placed below the belly button, 1 cm above the pubic bone (2) to detect the bladder (1).

Follow the directions in the app until the position is OK and the bladder has been found.


If the bladder is found, the position search light is green for 3 seconds. The app shows “bladder found.” If the bladder is not found within 1 min, the app shows “We couldn’t find your bladder”. Stop the position check in case the bladder is not full enough and repeat at least 45 minutes after last urination. The position check can be stopped at any time.

Mark your sensor location

Make a mark on your skin with a pen, 1 cm above the sensor. This will be used as a guide so you can place it correctly when the adhesive has been attached.

Note: The next time you use the sensor, you do not have to repeat the position check, but can put the sensor at this mark.

Remove ultrasound gel

Remove the gel from the skin and the back of the sensor.

Placing the adhesive patch

The adhesive patch needs to be replaced every day. Remove any excessive hair growth for an optimal bond of the adhesive patch to the skin.

Prepare adhesive patch

Take a patch and remove the upper foil layer.

Place adhesive patch

Place the patch on your sensor, aligning the opening on the patch with the black rectangle of the sensor.

Apply ultrasound gel

Apply gel to the black rectangular area on the back of your sensor.

Remove foil

Remove the last foil layer from the adhesive patch.

Make sure the sensor is attached to the correct position on the lower abdomen as determined during the position check. The position check can be repeated if desired.

Attach the Bladder Sensor

Attach the sensor to your abdomen, aligning it with the mark you made earlier. Placement done!

Step 4. Customize sensor settings

The sensor notifies the user when the bladder is full and nearly full and can provide a signal via the app, if desired.

- If the vibration is activated:
 The sensor will vibrate every 3 minutes until the bladder is emptied. After the first reminder the sensor will vibrate every 3 minutes with increased intensity.

- If app-notifications are activated:
 The app will provide a notification when the bladder is nearly full, and full.
- Battery too low:
 The app will provide a notification when the battery is below 15% and below 5%.

Your preferred type of notification can be adjusted in the settings of the app.

Select your type of reminder


- Just before you should go (when your bladder is 90% and 100% full)
- A while before you should go (when your bladder is 80% and 100% full)
- A good time before you should go (when your bladder is 70% and 100% full)

You can choose between a direct “your bladder is 100% full” or discrete “Time for a break!” notification. The preferred type of notification can be managed in the settings of the app.

4.3 WHILE USING THE BLADDER SENSOR


The app on the Smartphone displays the bladder filling status according to the last update.

Figure 1:
While using the sensor – Smartphone

The app on the Apple Watch displays the bladder filling status according to the last update.

Figure 2:
While using the Bladder Senso – Apple Watch

The overview in the app will indicate the level of bladder fullness with the following statuses:

LOW

Bladder is filled between 0% and 25%.


ABOVE 25%

The overview in the app will indicate the level of bladder fullness by showing a wave with a corresponding bladder filling percentage.

In addition, the overview in the app shows if the sensor is connected and how much power is left in the battery (see p. 25 and p. 26).

4.3.1 HISTORY

The sensor registers the voiding behavior of the user. The voiding diary can be accessed by selecting the 'history' page on the overview page in the app. The diary indicates the time the bladder is emptied and when the user has received a notification.

History


SUMMARY

This overview provides insight into the voiding pattern of the user.

- Indicates a voiding after a notification
- Manually specified that an accidental leakage occurred.

4.3.1.1 WHILE USING THE VOIDING DIARY

The voiding diary provides the information by:

WEEK

Overview presents your voiding behavior per week (one day per column).

Select "MONTH" on the voiding diary app to gain access to a detailed overview of the whole month.

MONTH

Overview presents your voiding behavior per month (one week per column).

en

Select “DETAILS” to gain access to a detailed overview of a specific day.

DAY

Overview presents your voiding behaviour per day. Indicating the time from the reminder until the bladder is emptied.

WARNING: New data for the voiding diary can only be provided to the app if the sensor is in proximity and connected to the Smartphone and/or Apple Watch. In-between different users, data can be cleared before reuse.

The distance between the sensor and the Smartphone and/or Apple Watch should not be more than ~10 meters in open air!

en

4.4 CHANGING SETTINGS

The app can be used to change the settings of how you want to be notified of a full bladder. Enter the notification settings menu by touching the settings (⚙) icon (see p. 25 and p. 26).

- Notifications on/off – This setting allows you to mute the sound of the text message you receive when the bladder is full.
- Notification type – This setting allows you to choose the desired type of notification; direct or discrete.
- Sensor vibration on/off – This setting allows you to mute the sound of the vibration, in total, of the device when the bladder is full. You can also adjust the intensity (low – medium – high).

WARNING: When the vibration is turned off, the user can only receive notification via the Smartphone and/or Apple Watch that is connected to the sensor. The connected Smartphone and/or Apple Watch must be within 10 meters of the user at all times!

4.5 HOW TO STOP USING THE BLADDER SENSOR

Removing the sensor after use

Press the skin while removing the sensor and the adhesive patch from the skin (replace the adhesive patch within 24 hours of use).

Clean the skin:

- Clean the abdomen with a tissue.

Clean the Bladder Sensor:

Step 1

- Remove the adhesive patch from the sensor. The adhesive patch can be thrown away as residual waste.

Step 2

For single patient, multiple use

- Clean the sensor with a soft, damp cloth after every use to remove any excess ultrasound gel, and wipe dry with a clean cloth.

- The following products may be used to clean the exterior surfaces of the sensor:
 - (Isopropyl) alcohol (70% solution in water)
 - Mild soap and water

Recharge the sensor for the next time.

For multiple patients or when the Bladder Sensor is clearly contaminated (e.g. urine or feces):

- Clean the sensor with a soft, damp cloth after every use to remove any excess ultrasound gel, and wipe dry with a clean cloth.
- The following products may be used to clean the exterior surfaces of the sensor:
- Clean by washing with lukewarm water and mild soap/detergent.
- Let airdry or wip dry with a paper towel and then clean the exterior surfaces with 70% alcohol (70% solution in water) and let again airdry before next usage.

Note: Use the USB charger / cable to re-charge the device (see chapter 1.1 section B and C)

WARNING: Make sure excess moisture stays away from openings. Never let the Bladder Sensor soak!

4.6 INDICATIONS THAT THE BLADDER SENSOR IS NOT WORKING PROPERLY

The Bladder Sensor is most probably not working properly if:

- The battery and position search lights are blinking red/orange.
- The battery light does not turn to green after charging for a couple of hours.
- The position light does not turn on after having given a position check command via the app.
- The information of the voiding diary is not available while connected to the sensor.

If:

- The app on the smartphone does not connect with the Bladder Sensor even when they are within 10 meters of each other, or the Apple Watch does not connect with the Bladder Sensor even when they are within 3 meters of each other; please restart the app on your Smartphone and/or Apple Watch.

In these cases, please stop usage and contact your local customer service: www.tena.com

Disposal:

- Expected lifetime and the warranty period of the TENA SmartCare Bladder Sensor is two years.
- The TENA SmartCare Bladder Sensor has been designed with recycling, reuse, and waste management in consideration. Novioscan will provide the necessary information to recycling facilities and government agencies as required by WEEE Directive.
- The adhesive patches and the ultrasound gel are materials which can be thrown away as residual waste (conform local regulations).

Malfunctioning and/or any serious incident:

- In case of malfunctioning and/or any serious incident that has occurred and is related to the Bladder Sensor, stop usage.

For more information in these cases, or in the case of other malfunctions, please report or contact the your local customer service via www.tena.com

Any serious incident should also be reported to the relevant competent authority of the EU Member State.

- In case of malfunctioning within the warranty term, your Bladder Sensor will be replaced with a well-functioning Bladder Sensor. If this is the case after the warranty term is expired, based on our assessment we shall suggest a mutually best-fitting solution.

6 CONTRAINDICATIONS

- The Bladder Sensor is not intended for pregnant or breastfeeding women.
- The Bladder Sensor is not intended for children with a BMI higher than the 85th percentile and adults with a BMI > 25 kg/m^2.
- The Bladder Sensor is not intended for monitoring bladder content in children < 30 ml, and in adults <100 ml (lower limits of sensitivity).
- The Bladder Sensor is not intended for notification purposes in children with a low night-time bladder capacity or bladder overactivity during the night (only use for monitoring purposes).
- The Bladder Sensor may only be used on intact skin.
- The Bladder Sensor is not intended for users in a wheelchair.

Read, understand and practice the precautionary and operating instructions.

- Only use the Bladder Sensor in combination with Novioscan's adhesive patches, USB charger and ultrasound gel.
- Minors must be aided by an adult during the setup and use of the Bladder Sensor.
- Handle the Bladder Sensor with care. Inappropriate handling of the Bladder Sensor may adversely affect its performance.
- Inspect the Bladder Sensor for cracks which may allow the ingress of conductive fluid before each use.
- DO NOT disassemble, modify or remodel the Bladder Sensor or its accessories. This may adversely affect its performance.
- DO NOT allow the Bladder Sensor to be heated above 45 °C during storage (113 °F).
- DO NOT allow the Bladder Sensor to be heated above 40 °C during operation (104 °F).
- DO NOT use in combination with X-ray or other ultrasound equipment.

- The USB charger (with USB charger cable) must always be able to be removed from a wall outlet / socket without any obstruction.
- The Bladder Sensor must not be immersed in water or other fluid.
- The Bladder Sensor may only be repaired, and the battery may only be replaced by qualified personnel (trained by Novioscan).
- The Bladder Sensor should not be used in an environment where Bluetooth® Low Energy is not allowed.
- Ensure that the environment in which the Bladder Sensor is used is not subject to any source of strong electromagnetic interference such as radio frequencies and magnetic fields. Potential electromagnetic or other interference could affect the Bladder Sensor or the other equipment. Try to minimize this interference: DO NOT use other equipment in conjunction with the Bladder Sensor.

7.1 ADDITIONAL PRECAUTIONARY INFORMATION FOR THE HOME ENVIRONMENT

- Do not lay the Bladder Sensor near the heater.
- Do not use the Bladder Sensor in the bathtub or under the shower.
- Do not soak the Bladder Sensor in water.
- Do not allow pets to chew or play with the Bladder Sensor or its accessories.
- Do not allow children to use the Bladder Sensor or its accessories for any other reasons than its intended use as described in section 4.1

TO PREVENT ELECTRICAL SHOCK, DO NOT OPEN THE SENSOR.

THERE ARE HIGH VOLTAGE COMPONENTS. LEAVE MAINTENANCE TO QUALIFIED STAFF ONLY.

en

7.2 CIRCUMSTANCES IN WHICH TO CONSULT A HEALTHCARE PROFESSIONAL

Consult a healthcare professional when use of the Bladder Sensor leads to:

- Skin irritation
- Painful skin during or shortly after use
- Damage to the skin
- Signs of infection such as redness of the skin, swelling, or blistering

en




8.1 TENA SMARTCARE BLADDER SENSOR

Ultrasound mode[1]**:** Ultrasound A-mode
Ultrasound frequency: 3.5 MHz
Sample frequency: 1.14 MHz
Measuring range: 1-173 mm
Battery: Lithium-ion battery DC 3.7V, 300mAh
Dimensions: ±75×61×14 mm
Weight: ±55 grams
Adhesive patch dimensions: 75×135 mm
Data memory: Flash drive
Data connection: Bluetooth® Low Energy 5.x to Smartphone app *(Inseparable part of the medical device: Bladder Sensor)*
Connector: USB Type C
Accessories: Power supply (USB charger & cable), ultrasound transmission gel, adhesive patches
Water resistance: IP65[2]

[1] *In accordance with IEC60601:2-37, the THERMAL INDICES and the MECHANICAL INDEX are below 1.0 for all device settings.*

[2] *Protected against access to a wire of 1.0 mm diameter, and no ingress of dust. Protected from water projected in jets against the enclosure from any direction.*

Environmental conditions during storage and use

	Ambient operating temperature: -10°C to +40°C
	Ambient temperature while charging: 0°C to +35°C
	Ambient temperature during storage: -10°C to +45°C
	Relative humidity for normal use and storage: 15%-75%, non-condensing
	Atmospheric pressure for normal use and storage: 700 to 1060 hPa

WARNING: If the TENA SmartCare Bladder Sensor is operating at a maximum environment temperature of 40°C and the 'Position Check' function is repeatedly activated, the sensor area might reach a maximum skin contact temperature of 43°C.

Remove the Bladder Sensor from the body if any significant temperature increase has been sensed and allow the Bladder Sensor to cool off for at least 30 minutes before attaching it to the body again. Strictly follow the Bladder Sensor precautionary information of chapter 7 for safe operation.

Essential Performance (IEC60601-2-37)[1]

The essential performance of the Bladder Sensor is:

- Free from the production of unintended or excessive ultrasound output.
- Free from the production of unintended or excessive TRANSDUCER ASSEMBLY surface temperature.

The Bladder Sensor is the actual applied part and is an internal battery powered device. The degree of shock protection is BF. The connection between the user and the device complies with permitted leakage currents and dielectric strengths of IEC/EN60601.

8.2 BLUETOOTH®

The sensor uses the 2.4 GHz (2.36 – 2.50 GHz) Bluetooth® module to connect with the Smartphone and/or Apple Watch. The Bluetooth® low energy (BLE) 5.3 is based on the nRF52840 Chip (modulation type). The maximum RF output power equals +8.0 dBm. The sensor should not be used in an environment.

8.3 USB CHARGER

REFERENCE	GTM46101-1005-USB
SIZE	41.0 x 71.0 x 31.5
OUTPUT VOLTAGE RANGE	5V
LOAD	Min: 0A, Max: 0.3A
RIPPLE (1) P-P (MAX)	150mV
REGULATION	Line: ±5%, Load: ±5%

8.3.1 INPUT


AC INPUT VOLTAGE RATING:	100 to 240VAC
AC INPUT VOLTAGE RANGE:	90 to 264VAC
MAX INPUT CURRENT	0.3A
INRUSH CURRENT	No damage. On power when switching on/off repeatedly under any condition

8.3.2 OUTPUT

MAX OUTPUT CURRENT	2A
MAX OUTPUT POWER RATING	10W
EFFICIENCY	EU CoC Ver. 5 Tier 2
INPUT POWER (NO LOAD)	< 0.1W < 0.1W, < 0.075W
HOLD-UP TIME	8mS Min @ 100VAC Full Load
OVER-VOLTAGE PROTECTION	Output shut down and auto restart 110% to 140% of Vnom

en

Label

Serial number (DM)

(01) Product code [14 digits]

(11) Production date [YYMMDD]

(21) Serial number [<15 digits>]

en

On the TENA SmartCare Bladder Sensor itself

SYMBOL	MEANING
	Manufacturer of the TENA SmartCare Bladder Sensor (filled).
MD	Medical device.
SN	Serial number.
	WEEE mark The device is electrical waste and should not be deposited among regular waste.
CE 1639	CE certified product: The product system conforms to the General Safety and Performance Requirements of the EU MDR 2017/745.
IP65	International Protection against contact and ingress of dust and against water jets.
	Applied Part, type BF (medium- or long-term contact with patient).

en

SYMBOL	MEANING
	Consult instructions for use (paper-based)
	The GS1 2D DataMatrix
N/A	Made in NL
N/A	Product name: TENA SmartCare Bladder Sensor
	TENA logo

Additional symbols on the packaging and in the Instructions For Use:

SYMBOL	MEANING
	Reference / article number
	Unique device identifier (UDI)
	Date of Manufacturing Format YYYY-MM-DD
	Radio Frequency Radiation Symbol




en

SYMBOL	MEANING
	Relative humidity requirements of environment during transport, storage, and use.
	Temperature requirements of the environment during transport, storage, and use.
	Atmospheric pressure requirements of the environment during transport, storage, and use.
	Warning
	Warning, electricity.
	Recycling mark Triman for FR-capsule
	Waste regulation-capsule Italy
	Waste regulation-plastic (due to ultrasound gel bottle) Italy
	**Distributor** Essity Hygiene and Health AB SE-405 03 Göteborg, Sweden Visiting address: Mölndals Bro 2, Mölndal www.essity.com
	Use-By Date (YYYY-MM-DD)


en

Shock-free, non-crushing packaging.

Environmental conditions during transport.
Bladder Sensor / Consumer Pack:

SYMBOL	MEANING
	Temperature during transport[1]: -10°C to +45°C
	Relative humidity during transport 15%-75%, non-condensing
	Atmospheric pressure: 700 to 1060 hPa

Bladder Sensor / Starter Kit:


SYMBOL	MEANING
	Temperature during transport[1]: +1°C to +40°C
	Relative humidity during transport 15%-75%, non-condensing
	Atmospheric pressure: 700 to 1060 hPa

[1] *If the device is transported outside operating temperature, it should be placed in an environment with the ambient temperature for 6 hours before used as intended. If the ambient temperature is reached earlier, the device can be used earlier.*


en

10.1 TRANSPORTATION LABELS

TENA SmartCare Bladder Sensor / Consumer Pack label:

TENA SmartCare Bladder Sensor / Starter Kit Label:

en

11.1 MANUFACTURER

Novioscan B.V.
Transistorweg 5
6534 AT Nijmegen
The Netherlands

en

⚠ WARNING

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

⚠ WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the sensor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

EMC tests have been passed under the following operation modes:

OPERATION MODE	DESCRIPTION
OM#01	EUT ON: Stand alone basis. Powered by batteries.
OM#02	EUT ON: Bluetooth® link established. Powered by batteries.

12.1 ELECTROMAGNETIC EMISSIONS

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS

The sensor is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the sensor is used in the following environment:

EMISSIONS TEST	COMP-LIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The sensor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

12.1 ELECTROMAGNETIC EMISSIONS

EMISSIONS TEST	COMP-LIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Class B	The sensor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	Not applicable. Rated power is less than 75W
Voltage Fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

12.2 ELECTROMAGNETIC IMMUNITY

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY

The sensor is intended for use in the electromagnetic environment specified below. The customer or the user of the sensor should assure that the device is used in such an environment:

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMP-LIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be made of wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s)	±1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_T for 0,5 and 1 cycle 70% U_T for 25/30 cycles 0% U_T for 250/300 cycles	0% U_T for 0,5 and 1 cycle 70% U_T for 25/30 cycles 0% U_T for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the sensor requires continued operation during power mains interruptions, it is recommended that the sensor be powered from an uninterruptible power supply or a battery.

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMP-LIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT GUIDANCE
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T IS THE A.C. MAINS VOLTAGE PRIOR TO APPLICATION OF THE TEST LEVEL

RF Common mode/ Conducted Susceptibility IEC 61000-4-6	3 Vrms 6 V in ISM and amateur radio bands 150 kHz to 80 MHz	3 Vrms 6 V in ISM and amateur radio bands	

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMP-LIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT GUIDANCE
Radiated RF Electromagnetic Field IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m	
Test specifications for enclosure port immunity to proximity magnetic fields IEC 61000-4-39 (Table 11, IEC 60601-1-2:2014+A1: 2020)	See next table below	See next table below	See next table below

TEST FREQUENCY	MODULATION	IMMUNITY TEST LEVEL (A/M)
30 kHz a)	CW	8
134,2 kHz	Pulse modulation b) 2,1 kHz	65 c)
13,56 MHz	Pulse modulation b) 50 kHz	7,5 c)

a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.
b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
c) r.m.s., before modulation is applied.

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT GUIDANCE
Proximity fields from RF wireless communication systems IEC 61000-4-3	See next table	See next table	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the sensor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Interference may occur in the proximity of equipment with the following symbol: 

TEST FREQUENCY (MHZ)	BAND (MHZ) a)	SERVICE a)	MODULATION b)	MAXIMUM POWER (W)	DISTANCE (M)	IMMUNITY TEST LEVEL (V/M)
385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704-787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
870						
930						
1720				2	0,3	28

TEST FREQUENCY (MHZ)	BAND (MHZ) a)	SERVICE a)	MODULATION b)	MAXIMUM POWER (W)	DISTANCE (M)	IMMUNITY TEST LEVEL (V/M)
1 845 – 1 970	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz			
2 450	2 400 – 2 570	Bluetooth®, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
5 240 – 5 500 – 5 785	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The regulations applicable to the TENA SmartCare Bladder Sensor are:

- EU MDR 2017/745
 Medical Device Regulation (MDR)
- 2014/53/EU
 Radio Equipment Directive (RED)
- 94/62/EC
 Packaging and packaging waste
- 2011/65/EU
 Restriction of hazardous substances (RoHS)
- 2012/19/EU
 Waste from electrical and electronic equipment (WEEE)

Using the list above as a guide and after consultation with DEKRA's test house, the standards in the sections below have been found applicable to the sensor.

In compliance with EU MDR 2017/745:

NORM	TITLE
ISO 15223-1:2021	Medical devices — Symbols to be used with medical device labels, labeling and information to be supplied.
NEN-EN-ISO 13485:2016	Medical devices — Quality management systems — Requirements for regulatory purposes.
NEN-EN-ISO 14155:2020	Clinical investigation of medical devices for human subjects — Good clinical practice.
NEN-EN-ISO 14971:2019	Medical devices — Application of risk management to medical devices.
IEC 60601-1:2005 + A1:2012 + A2:2022	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.

NORM	TITLE
IEC 60601-1-11:2015 + A1:2020	Medical electrical equipment — Part 1-11: General requirements for basic safety and essential performance — Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment).
IEC 60601-2-37:2007 + A1:2015	Medical electrical equipment — Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment.
EN 60601-1-2:2014 + A1:2020	Medical electrical equipment — Part 1-2: General requirements for basic safety and essential performance. Collateral standard: Electromagnetic compatibility.
IEC 62304: 2006+A1:2015	Medical device software — Software life cycle processes.

NORM	TITLE
IEC 62366-1:2015 + A1:2020	Medical devices — Application of usability engineering to medical devices.
ISO 20417:2021	Medical devices — Information to be supplied by the manufacturer.
IEC 62133:2017	Secondary cells and batteries containing alkaline or other non-acid electrolytes — Safety requirements for portable sealed secondary lithium cells, and for batteries made from them, for use in portable applications — Part 2: Lithium systems.
ISO 10993-1:2018	Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management.
ISO 10993-5: 2009	Biological evaluation of medical devices — Part 5: Tests for in vitro cytotoxicity.
ISO 10993-10:2021	Biological evaluation of medical devices — Part 10: Tests for skin sensitization

In compliance with RED 2014/53/EU:

NORM	TITLE
EN 301 489-1 V2.1.1 (2017) /V1.9.2	ElectroMagnetic Compatibility (EMC) standards for radio equipment and services; Part 1: Common technical requirements; Harmonised Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU and the essential requirements of article 6 of Directive 2014/30/EU.
EN 301 489-17 V3.2.4 (2020-09)	ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems; Harmonised Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU.
EN300 328 v2.2.2 (2019-07)	Wideband transmission systems; Data transmission equipment operating in the 2.4 GHz ISM band and using wide band modulation techniques; Harmonised Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU.

NORM	TITLE
IEC62479:2010	Assessment of the compliance of low-power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300GHz).

In compliance with Packaging and packaging waste 94/62/EC:

NORM	TITLE
EN 2248:1993	Packaging - Complete, filled transport packages - Vertical impact test by dropping.
ISTA6	Ships in Own Container (SIOC) for Amazon.com Distribution System Shipment.

In compliance with Restriction of hazardous substances (RoHS) 2011/65/EU:

NORM	TITLE
EN 50581:2012	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances.

The requirements of the Directive 'Waste from electrical and electronic equipment (WEEE) 2012/19/EU' have been followed.

- The proper disposal of WEEE products will conserve natural resources (through reuse and recycling) and protect environmental and human health through proper disposal.
- WEEE-labelled products must be kept separate and returned to a designated collection point for dismantling, recycling and/ or reuse.
- The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local government, your commercial, institutional and/ or municipal waste disposal service or the location where you purchased the product.

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