

NOT ALL ARMOUR IS MADE OF STEEL



JOBST

DEFY GRAVITY

JOBST® FarrowWrap®

Simple decongestion therapy through adjustable
short stretch wrap solutions


FITTED FOR LIFE

At Essity, we focus on integrated therapy solutions and so caring for lymphedema means supporting patients every step of the way. From our range of daytime and nighttime custom and ready-to-wear products, to our wide selection of bandages & wraps, our portfolio is highly adaptable which makes it easier for patients to take an active role in finding the therapies that are best suited for them.

DECONGESTION



► Reducing the edema to a manageable size and improving skin health



MAINTENANCE



► Keeping edema stable over time by adapting therapy to new needs



JOBST® FARROWWRAP®

Short stretch compression wraps

- Easy-to-use adjustable bands
- Easy to don
- For patients requiring independence



JOBST® COMPRI2

2-layer short stretch compression systems

- Sustained compression for up to 7 days*
- Simple application with visual compression indicators



COMPRILAN®

Short stretch compression bandages

- Simple and cost-effective solution
- Ideal when frequent applications are required

DAYTIME



JOBST® ELVAREX®

Custom-made compression garments

- Gold standard in lymphedema management
- Highly individualizable products for upper and lower extremities
- Unique features for better comfort and fit (Functional Zones and SoftFit™ Technology)

*Caution: These products are made with natural rubber latex.

NIGHTTIME



JOBST® JOVIPAK®

Custom-made compression garments

- Promotion of lymph drainage
- Softening of fibrosis
- Increased patient independence

*Targeted nominal compression for an average ankle size.

ENABLING LYMPHEDEMA CARE EVERY DAY WITH JOBST® FARROWWRAP®

We consulted patients with lymphedema and have seen how difficult therapy can be. Patients struggle to don garments and wear them comfortably throughout the day. They are also challenged by time-consuming therapies such as bandaging.

WHO BENEFITS FROM JOBST® FARROWWRAP®?

- Patients who require independence
- Patients who find compression bandaging difficult or that don't have regular access to a professional who can assist
- Patients with swelling that fluctuates




HOW DOES JOBST® FARROWWRAP® WORK?

- It provides steady and comfortable resting compression, but higher working compression when the patient is active.
- Helps prevent limb from swelling during the day or night.



JOBST® FARROWWRAP® PRODUCT SELECTION

1. CHOOSE COMPRESSION LEVEL AND FABRIC

	Fabric*	For who?	Features	Variants
Mild to Moderate 20-30 mmHg	LITE	Neuropathic and pediatric patients or in palliative care	- Double laminated fabric - Conforms to limb shape	 RTW, TTF, CM
	CLASSIC	- Original fabric - Suitable for stubborn edema with skin folds and irregular limb shapes	Durable, rigid fabric with a comfortable feel	 RTW (except arm), CM
	STRONG	- Suitable for minor skin folds and distortions - Less rigid than CLASSIC with a soft inner layer	- Contoured bands conforms to limb - Highly breathable - Double laminated fabric	 RTW (except arm), TTF, CM
	BASIC	Cost-effective fabric with quality compression	Double-sided VELCRO® so you can change band length and trim as limb reduces	 RTW
	4000	Easy to apply with just four bands	Inner sleeve assists donning	 RTW

*All fabrics are latex free

RTW: Ready-To-Wear • CM: Custom-Made • TTF: Trim-To-Fit

2. MEASURE TO CHOOSE STYLE

Product	Description
READY-TO-WEAR (RTW)	• Range of pre-made sizes • More cost-effective than custom garments
TRIM-TO-FIT (TTF)*	• Can be trimmed to adapt as extremity reduces • Can fit most patients in a single size • Upper and lower extremities
CUSTOM-MADE (CM)*	• Customized garments for patients who cannot properly fit in a ready-made garment

* An order form must be completed for Trim-to-Fit and Custom-Made garments

3. CHOOSE LINER AND ACCESSORIES**

Product	Description
JOBST® FarrowHybrid® ADI Foot Compression Liner	• 20-30 mmHg • Alternative to using a footpiece
Flesh Tone Liner	• No compression • Adjusts to calf up to 58 cm circumference
TG® Soft Padded lining	• No compression • Tubular shape • Provides extra comfort and padding • Upper and lower extremities
JOBST® Farrow LITE Gauntlet	• 20-30 mmHg • Ambidextrous hand gauntlet

**A liner is required for all upper and lower limb pieces with the exception of the RTW armpiece

ARMPIECE



JOBST® FarrowWrap® LITE

Description	Colour	Length	XS	S	M	L	XL
RTW Armpiece Right	Tan	Short		7665012	7665014	7665016	
RTW Armpiece Left	Tan	Short		7665011	7665013	7665015	
RTW Armpiece Right	Tan	Regular		7665006	7665008	7665010	
RTW Armpiece Left	Tan	Regular		7665005	7665007	7665009	
RTW Armpiece Right	Tan	Long		7665000	7665002	7665004	
RTW Armpiece Left	Tan	Long		7665035	7665001	7665003	
RTW Armpiece Right	Black	Short		7665030	7665032	7665034	
RTW Armpiece Left	Black	Short		7665029	7665031	7665033	
RTW Armpiece Right	Black	Regular		7665024	7665026	7665028	
RTW Armpiece Left	Black	Regular		7665023	7665025	7665027	
RTW Armpiece Right	Black	Long		7665018	7665020	7665022	
RTW Armpiece Left	Black	Long		7665017	7665019	7665021	
Ambidextrous Hand Gauntlet with Foam Squares	Tan	Regular	7665300	7665301	7665302	7665303	
TTF Armpiece Right (includes Velcro)	Tan	Regular			7612002		7612003
TTF Armpiece Left (includes Velcro)	Tan	Regular			7612000		7612001
TTF Armpiece Extra Velcro	7612004						

SIZING CHART

Measurement (cm)	S	M	L
G - Axilla circumference	22-31	29-39	32-45
E - Elbow Crease circumference	20-27	25-34	30-40
C - Wrist Crease circumference	14-18	16-21	19-25
C-G Short*	40-43		
C-G Regular*	44-47		
C-G Long*	48-51		

*C-G: Length from the wrist to the axilla



LITE RTW



CLASSIC Custom



LITE Custom
STRONG Custom
LITE TTF (adjustable)

THIGHPIECE



JOBST® FarrowWrap® LITE

Description	Colour	Length	XS	S	M	L	XL
RTW Thighpiece	Tan	Short	7665505	7665506	7665507	7665508	7665509
RTW Thighpiece	Tan	Regular	7665500	7665501	7665502	7665503	7665504
RTW Thighpiece	Tan	Long	7665510	7665511	7665512	7665513	7665514

JOBST® FarrowWrap® CLASSIC

Description	Colour	Length	XS	S	M	L	XL
RTW Thighpiece	Tan	Short	7665806	7665807	7665808	7665809	7665810
RTW Thighpiece	Tan	Regular	7665800	7665801	7665802	7665803	7665804
RTW Thighpiece	Tan	Long	7665812	7665813	7665814	7665815	7665816

JOBST® FarrowWrap® STRONG

Description	Colour	Length	XS	S	M	L	XL
RTW Thighpiece	Tan	Short	7664905	7664906	7664907	7664908	7664909
RTW Thighpiece	Tan	Regular	7664900	7664901	7664902	7664903	7664904
RTW Thighpiece	Tan	Long	7664910	7664911	7664912	7664913	766491

SIZING CHART

Measurement (cm)	XS	S	M	L	XL
Groin circumference	60-70	67-75	70-80	75-90	85-100
Above Knee circumference	47-52	53-59	60-66	67-73	74-80
E-G Short*	16-20				
E-G Regular*	21-26				
E-G Long*	27-31				

*E-G: Medial length from just above patella to the groin. Follow contour of any skin folds/lobules.
Most thighpieces can accommodate a 20% reduction of the circumference.



LITE
STRONG



CLASSIC

LEGPiece



JOBST® FarrowWrap® LITE

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7665400	7665401	7665402	7665403	7665404
RTW Legpiece	Tan	Long	7665405	7665406	7665407	7665408	7665409
TTF Legpiece (includes Velcro)					7612300		
TTF Legpiece Extra Velcro	7612301						
TTF Legpiece Extra Band	Tan	Regular			7612302		

JOBST® FarrowWrap® CLASSIC

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7665700	7665701	7665702	7665703	7665704
RTW Legpiece	Tan	Long	7665705	7665706	7665707	7665708	7665709

JOBST® FarrowWrap® STRONG

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7664800	7664801	7664802	7664803	7664804
RTW Legpiece	Tan	Long	7664805	7664806	7664807	7664808	7664809
TTF Legpiece (includes Velcro)	Tan	Regular			7611800		7611801
TTF Legpiece Extra Velcro	7611802						
TTF Legpiece Extra Band	Tan	Regular			7611803		7611804

JOBST® FarrowWrap® BASIC*

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece (includes Velcro)	Tan	Regular	7666000	7666001	7666002	7666003	7666004
RTW Legpiece (includes Velcro)	Tan	Long	7666005	7666006	7666007	7666008	7666009
RTW Legpiece Extra Velcro	7666010						

*Comes with a JOBST® FarrowHybrid® ADI Wide Foot Compression Liner

LEGPiece

LEGPiece SIZING CHART (except JOBST® FarrowWrap® 4000)

Measurement (cm)	XS	S	M	L	XL
Mid-Calf circumference	36-43	42-50	48-58	53-63	58-68
Ankle circumference	21-25	25-30	30-36	36-42	42-50
A-D Regular*	33-37	35-39	37-41	39-43	39-43
B-D Regular**	28	30	31	33	33
A-D Long*	38-41	40-43	42-45	44-47	44-47
B-D Long**	33	35	36	38	38

Circumferences measurements must be within ± 1.0 to 1.5 cm

*A-D: Floor length up to 2 fingers below knee fold

**B-D: Actual height of the garment

JOBST® FarrowWrap® 4000*

Description	Colour	Length	XS	S	M	L
Legpiece	Tan	Regular	7666100	7666101	7666102	7666103
Legpiece	Tan	Long	7666104	7666105	7666106	7666107
Legpiece	Black	Regular	7666108	7666109	7666110	7666111
Legpiece	Black	Long	7666112	7666113	7666114	7666115

*Comes with a JOBST® FarrowHybrid® ADI Wide Foot Compression Liner

SIZING CHART (JOBST® FarrowWrap® 4000)

Measurement (cm)	XS	S	M	L
Mid-Calf circumference	25-38	32-48	37-56	43-67
Ankle circumference	17-23	19-28	25-37	29-43
B-D Regular*	32			
B-D Long*	36			

*B-D: Actual height of the garment



LITE
STRONG
BASIC



CLASSIC



4000

FOOTPIECE



JOBST® FarrowWrap® LITE

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7665205	7665206	7665207	7665208	7665209
RTW Legpiece	Tan	Long	7665200	7665201	7665202	7665203	7665204
TTF Footpiece (includes Velcro)	Tan	Regular			7612100		
TTF Footpiece Extra Velcro	7667500						

JOBST® FarrowWrap® CLASSIC

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7665605	7665606	7665607	7665608	7665609
RTW Legpiece	Tan	Long	7665600	7665601	7665602	7665603	7665604

JOBST® FarrowWrap® STRONG

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7664604	7664605	7664606	7664607	7664608
RTW Legpiece	Tan	Long	7664609	7664600	7664601	7664602	7664603
TTF Footpiece (includes Velcro)	Tan	Regular			7611600		
TTF Footpiece Extra Velcro	7611601						

JOBST® FarrowWrap® BASIC*

Description	Colour	Length	XS	S	M	L	XL
RTW Legpiece	Tan	Regular	7665905	7665906	7665907	7665908	7665909
RTW Legpiece	Tan	Long	7665900	7665901	7665902	7665903	7665904

*Comes with a JOBST® FarrowHybrid® ADI Wide Foot Compression Liner



BASIC

SIZING CHART

Measurement (cm)	XS	S	M	L	XL
Mid-foot circumference	22-24	25-27	28-30	31-34	35-40
Regular*	16-17	17-18	19-20	20-21	22-23
Long**	18-19	19-20	21-22	22-23	24-25

*Circumferences measurements must be within ± 0.5 cm

**Length from back of the heel to base of pinky toe



CLASSIC



LITE
STRONG

ACCESSORIES

JOBST® FarrowHybrid®

Description	Colour	S	M	L	XL
ADI Regular	Tan	7666400	7666401	7666402	7666403
ADI Large	Tan	7666412	7666413	7666414	7666415
ADI Large	Black	7666416	7666417	7666418	7666419

Calf Circumference must be less than 60 cm

SIZING CHART

Measurement (cm)	Regular	Large
Mid-foot circumference	20-27	28-35

	S	M	L	XL
Men Shoe Size	3-6	7-9	10-12	13-15
Women Shoe Size	4-7	8-10	11-13	14-16

*Do not apply the JOBST® FarrowHybrid® if calf circumference exceeds 60cm.
Smooth out any wrinkles in the garment.

Flesh Tone Liners

All Flesh Tone and TG® Soft liners are sold as pairs. Free Knee High liner issued with any legpiece or footpiece & legpiece combination.

Description	Colour	S	M	L
Knee High	Tan	7666621		7666622
Arm	Tan	7666620		

TG® Soft Liners

Description	Colour	S	M	L
Knee High	Tan	7666900	7666902	7666904
Thigh High	Tan	7667000	7667002	7667004
Arm	Tan	7667006	7667008	7667010

SIZING CHART

Measurement (cm)	S	M	L
TG® Soft			
Arm, calf or groin circumference	<40	41-70	71-120
Flesh Tone Liners			
Knee High	19-48 cm		48-70 cm
Arm	19-48 cm		48-70 cm

ACCESSORIES

JOBST® FarrowWrap® LITE Gauntlet

Description	XS	S	M	L
Gauntlet	7665300	7665301	7665302	7665303



LITE Gauntlet

SIZING CHART

Measurement (cm)	XS	S	M	L
Base of fingers circumference	16-18	18-20	20-22	22-24
Fold of thumb circumference	17-19	19-21	21-23	23-26
Wrist circumference	15-18	17-20	19-22	22-24

WASH AND CARE GUIDE

Description	CLASSIC	LITE	STRONG	BASIC	Hybrid ADI	Flesh Tone & TG® Soft	GarmentGrip™
Hand Wash	X						
Machine Wash Warm					X	X	
Machine Wash Cold		X	X	X			X
Lay Flat to Dry	X						
Tumble Dry, Low Heat					X	X	
Tumble Dry, No Heat		X	X	X			
Do Not Bleach	X	X	X	X	X	X	X
Do Not Iron	X	X	X	X	X	X	X
Do Not Dry Clean	X	X	X	X	X	X	X