

Ready-to-Wear

**Armsleeves and Gauntlets** 

### JOBST® Bella™ Strong

#### **Comfortable Compression**

#### Improved comfort, better compliance, and enhanced quality of life

Patients are more likely to be compliant when they feel comfortable wearing a compression garment. **JOBST® Bella™ Strong Ready-to-Wear** products – worn together or separately – effectively control for moderate lymphedema and edema.

#### JOBST® Bella™Strong Armsleeve

- Comfort Flat wristband reduces rolling and improves durability
- · Available with or without silicone band
- Easy to don, armsleeve comes with donning device
- · Available with knitted top band

#### JOBST® Bella™ Strong Gauntlet

- Avoids webbing at thumb area
- Extended wrist avoids pinching at the wrist
- · Comfort Flat wrist seam reduces rolling

### JOBST® Bella™ Lite

#### The Entry Point to JOBST® Quality

The **JOBST**® **Bella**™ **Lite** Ready-to-Wear Armsleeve and Gauntlet, together or separately, help maintain the limb size of lymphedema patients. These comfortable garments provide gradient compression to assist in patient management following trauma or surgery.

- Soft and air permeable, lightweight, durable knit construction
- Reduced compression at the wrist provides an effective and comfortable fit when the Armsleeve and Gauntlet are worn together
- Easy to put on, comfortable to wear
- Armsleeve available with silicone dot band



#### INDICATIONS JOBST® BELLA™ STRONG

Armsleeve Compression garments can also be used for: Lymphedema and other edema

Compression garments can also be used for:

Phlebitis • Post-thrombotic syndrome • Vascular malformations

#### INDICATIONS JOBST® BELLA™ LITE

Armsleeve Compression garments can also be used for: Lymphedema and other edema

Compression garments can also be used for:

Phlebitis • Post-thrombotic syndrome • Vascular malformations • Post-liposuction

#### **CONTRAINDICATIONS**

Arterial insufficiency • Untreated septic phlebitis
Conditions where increased venous and/or lymphatic return is not desired

## **COMFORT** for you

#### **Comfort and Performance**

JOBST® Bella™ Strong provides patients with the latest development in fiber science. The skin-friendly multifiber nylon combines the clinically effective management of lymphedema with a natural feeling all day long, day after day.



JOBST® Advanced comfort

skin

perspiration

evaporation

Quickly wicks moisture

Keeps you cool and comfortable

Feels soft to the skin

Helps to regulate relative

humidity, feeling cooler in warm conditions and warmer in cool

conditions

In addition, the open knit structure allows for an improved air permeability that keeps patients' skin drier than the competition.



Bella™ Strong: 20-30 mmHd



Bella<sup>™</sup> Strong:  $30-40\,\text{mmHg}$ 



Competitor I: 30-40 mmHg



Competitor II: 30-40 mmHg

## **INNOVATION** for you


#### Innovative sizing system - Covers an unrivalled range of patient dimensions

Considering over 50 years of scientific research and practical experience as well as the data from over 36,000 lymphedema patients for the armsleeve and over 33,000 for the gauntlet, JOBST ® has developed an innovative measuring-size-system.

Catering to all relevant armsleeve measuring points – wrist, mid-lower and mid-upper arm – this thorough approach, accepts the largest range of patient dimensions. **JOBST® Bella™ Strong** armsleeve sizing begins at the forearm, as data shows that we can accommodate more patients; this allows for multiple size options in the wrist and mid-upper arm for greater coverage, recognizing that every patient is unique.


### JOBST® Bella™ Strong Armsleeve: Ten-size system provides over 80% fit coverage

Competitors' wrist centric size determination limits fit coverage.



### JOBST® Bella™ Strong Gauntlet: Over 90% fit coverage

With the eight-size system the **JOBST®Bella™ Strong** gauntlet reaches an outstanding range of fit coverage.



<sup>\*</sup> Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.

#### Donning made easy

It requires less force to put on the **JOBST® Bella™ Strong** Ready-to-Wear Armsleeve than those of the competition.

#### 

#### **Donning Aids**

**DYCEM® NON SLIP MATERIAL** non-slip surface for easier donning and improves grip.

JOBST® DONNING GLOVES\* create friction, provide grip and help to distribute the fabric evenly.



<sup>1</sup> Data on file

<sup>\*</sup>Caution: This product contains natural rubber latex which may cause allergic reactions.

## **SIZING** for you

#### JOBST® Bella™ Strong Armsleeve

Size	Wrist (cC)	Mid forearm (cD)	Mid upper (cG)
1	14-18 cm	17-21 cm	22-31 cm
2	14-18 cm	17-21 cm	31-41 cm
3	14-18 cm	21-25 cm	25-34 cm
4	16-19 cm	21-25 cm	28-38 cm
5	17 - 20 cm	21-25 cm	34-43 cm
_			
6	14-18 cm	25-29 cm	29-38 cm
7	17 - 20 cm	25-29 cm	33-43 cm
8	18-21 cm	25-29 cm	38-48 cm
9	16-19 cm	29-33 cm	32-42 cm
10	19-22 cm	29-33 cm	36-46 cm

#### **Sizing instructions**

Finding the right size for each patient is easy since stepwise exclusion leads to the perfect solution:

- 1. Measure wrist circumference.
- 2. Measure mid forearm circumference.
- 3. Measure circumference of the mid upper arm.

If you are left with two sizes to choose from, we recommend selecting the larger size for ease of donning, overall comfort, and compliance. Either size selection will deliver the desired therapy.

#### JOBST® Bella™ Strong Gauntlet

	Size	Palm at base of fingers (cA)	Wrist (cC)
$\Omega$	1	15 - 17 cm	14-18 cm
UMM	2	17 - 19 cm	14-18 cm
CA /	3	17 - 19 cm	18 - 21 cm
	4	19 - 21 cm	14-18 cm
	5	19 - 21 cm	18 - 21 cm
cC (	6	21-23 cm	16-19 cm
1 \	7	21-23 cm	19-22 cm
1 1	8	23-26 cm	19-22 cm


#### **Sizing instructions**

- 1. Measure palm circumference. This leaves 1 or 2 choices.
- 2. Measure wrist circumference.
- 3. Select product.




#### JOBST® Bella™ Lite Armsleeve

Size	Wrist	Mid-Lower	Mid-Upper
S	13.5 cm - 19.5 cm	17.8 cm - 23.2 cm	21.5 cm - 39.5 cm
М	15.5 cm - 21.5 cm	23.3 cm-28.7 cm	27 cm - 45 cm
L	17.5 cm - 23.5 cm	28.8 cm-34.3 cm	32.5 cm - 50.5 cm



All sizes are available in 2 lengths:

Regular (up to 43 cm); Long (greater than 43 cm)



#### JOBST® Bella™ Lite Gauntlet

	Size	Palm at base of fingers (cA)	Wrist (cC)
M a	s	15 cm - 19 cm	14 cm - 20 cm
cA	М	19 cm -22 cm	16 cm-22 cm
cC (	L	22 cm - 26cm	18 cm-24 cm

## **ORDERING** for you

#### JOBST® Bella™ Strong Armsleeve

**Compression Classes** 

Style

JOBST Bella™ Strong Armsleeve is available in ten ready-to-wear sizes.

#### Colour



30-40 mmHg

15-20 mmHg

20-30 mmHg

	15-20 mmHg* NATURAL	20-30 mmHg* NATURAL	30-40 mmHg* NATURAL
Size	Silicone Band	Silicone Band	Silicone Band
1	7272101	7273201	7277001
2	7272102	7273202	7277002
3	7272103	7273203	7277003
4	7272104	7273204	7277004
5	7272105	7273205	7277005
6	7272106	7273206	7277006
7	7272107	7273207	7277007
8	7272108	7273208	7277008
9	7272109	7273209	7277009
10	7272110	7273210	7277010

Please refer to your Price List for a complete list of catalogue numbers for ordering.

### JOBST® Bella™ Strong Gauntlet

**Compression Classes** 

15-20 mmHg

20-30 mmHg

Style



JOBST Bella™ Strong Gauntlet is available in eight ready-to-wear sizes.

Colour



	15-20 mmHg**	20-30 mmHg**
Size	Natural	Natural
1	7277600	7369410
2	7277601	7369411
3	7277602	7369412
4	7277603	7369413
5	7277604	7369414
6	7277605	7369415
7	7277606	7369416
8	7277607	7369417



<sup>\*</sup> Mean compression for average forearm size.

<sup>\*\*</sup> Mean compression for an average palm size.

# ORDERING for you

#### JOBST® Bella™ Lite Ready-to-Wear Armsleeve



#### Compression Classes

15-20 mmHg

20-30 mmHg

#### Style



JOBST Bella™ Lite Armsleeve is available in three ready-towear sizes and two lengths.

#### Colours





Size	Length	Beige	Black	mmHg
C	Regular	7714910		15 30 mm.
S	Long	7714901		15-20 mmHg
М	Regular	7715910		15 - 20 mm 4 a
IVI	Long	7715901		15-20 mmHg
L	Regular	7719110		15-20 mmHg
	Long	7719101		15-20 IIIIIIH
S/Silicone	Regular	7714911		15-20 mmHg
dot band	Long	7714902		13-20 111111119
M/Silicone	Regular	7715911		15-20 mmHg
dot band	Long	7715902		13-20 111111119
L/Silicone	Regular	7719111		15-20 mmHg
dot band	Long	7719102		15 20 11111119
S	Regular	7860810		20-30 mmHg
3	Long	7860804	7860809	20 30 11111119
M	Regular	7861310		20-30 mmHg
IVI	Long	7867304	7861309	20 30 11111119
L	Regular	7861410		20-30 mmHg
L	Long	7861404	7861409	20 30 11111119
S/Silicone	Regular	7860811		20-30 mmHg
dot band	Long	7860805		20 30 11111119
M/Silicone	Regular	7861311		20-30 mmHg
dot band	Long	7861305		20 30 IIIIIIII
L/Silicone	Regular	78661411		20-30 mmHg
dot band	Long	7861405		20 30 mining

#### JOBST® Bella™ Lite Ready-to-Wear Gauntlet



Size	Beige	Beige Black	
S	7719201		15-20 mmHg
М	7719901		15-20 mmHg
L	7727201		15-20 mmHg
S	7860904	7860907	20-30 mmHg
М	7863604	7863607	20-30 mmHg
L	7863704	7863707	20-30 mmHg

sizes.

Compression Classes

15-20 mmHg 20-30 mmHg Style



Colours JOBST Bella™ Lite

Gauntlet is available

in three ready-to-wear





black

## JOBST® Custom Flat-Knit Garments for Upper Extremities

When Ready-to-Wear sizing is not a perfect fit, customized garments are available through Elvarex®

#### JOBST® Elvarex® Custom Armsleeve / Glove / Gauntlet

CCL	Compression mmHg	Beige	Black	Navy	Caramel	Cranberry	Cherry	Bronze	Hazelnut
1	15 - 21 mmHg		•	•		•			•
2	23-32 mmHg		•	•		•			•
2 Forte	23-32 mmHg	•	•	•		•	•	•	

#### JOBST® Elvarex® Plus Custom Armsleeve / Glove / Gauntlet

CCL	Compression mmHg	Beige	Black	Navy	Caramel	Cranberry	Cherry	Red	Bronze	Hazelnut
1	15 - 21 mmHg			•		•				
2	23-32 mmHg		•	•		•		•		
2 Forte	23-32 mmHg		•	•		•		•		

#### JOBST® Elvarex® Soft Custom Armsleeve

CCL	Compression mmHg	Beige	Black	Navy	Honey	Cherry	Cocoa	Grey	Cranberry
1	15 - 21 mmHg		•	•		•	•	•	•
2	23-32 mmHg		•				•	•	

#### JOBST® Elvarex® Soft Seamless Custom Glove / Gauntlet

CCL	Compression mmHg	Beige	Black	Navy	Honey	Cocoa	Cherry	Grey	Cranberry
1	15 - 21 mmHg		•	•		•	•		•
2	23-32 mmHg		•	•		•		•	

Please refer to the product label and/or package insert for full instructions on the safe use of these products. All compressions shown are the mean compression for an average forearm size (armsleeve) or mean compression for an average palm size (glove / gauntlet).



