

JOBST®


The Complete Solution for Lymphedema

A Comprehensive Range of **Compression Products** designed with the **patient** and **therapist** in mind.

The Complete Solution for Lymphedema

We aim to support you and your patients along the full continuum of care with innovative products that provide therapeutic outcomes through improved patient compliance.

Why JOBST®?	3	JOBST® Custom Seamless Soft	
Cutimed® Siltec®		JOBST® Relax	10
Cutimed® Siltec® Sorbact®		JOBST® Bellavar®	
Cutimed® Sorbion®		JOBST® Custom Seamed	11
Cutimed® Sorbion® Sorbact®	4	JOBST® Bella™ Lite	
Compression Bandages		JOBST® Bella™ Strong	12
JOBST® Compri2	5	JOBST® SoftFit	
JOBST® FarrowWrap®	6	Band Variations	13
JOBST® Elvarex®	7	Functional Zones for JOBST® Elvarex®	
JOBST® Elvarex Plus®		and JOBST® Elvarex® Soft	14
JOBST® Elvarex® Soft	8	JOBST® Donning Aids & Accessories	15
JOBST® Confidence		JOBST® Educational Resources / Custom e-Shop	16
JOBST® Elvarex® Soft Seamless	9		



Why JOBST®?

Trust, confidence and history



- Conrad Jobst, mechanical engineer and industry pioneer, dedicated his life and expertise to inventing graduated compression garments in 1950.
- JOBST® is a brand of medical compression legwear for the prevention, treatment and management of venous and lymphatic conditions.
- After his passing, Conrad's wife Carolyn, assumed the presidency of the company for the next 27 years. Inspired by her late husband's legacy, Carolyn grew the business into a global organization and extended sales to virtually every continent as the JOBST® reputation for quality grew.

Innovation



- JOBST® incorporates customer and healthcare professional involvement in the development of products. From focus groups to wash-and-wear tests, all garments go through a rigorous design process that ensures they provide the highest efficacy, safety, and satisfaction for our customers.
- Innovation is at the heart of our product excellence. Product development efforts are directed towards delivering the highest quality product with the highest assurance of medical efficacy.

Partnerships



- JOBST® is an exclusive sponsor of the JOBST® Research Grant in Venous and Lymphatic Diseases in association with the American Venous Forum Foundation.
- JOBST® supports the National Lymphedema Network and Lymphatic Education and Research Network.
- Key JOBST® stakeholders are part of the Lymphedema Advocacy Group. Their mission is to advocate for improved insurance coverage for the diagnosis and treatment of the disease.



Education

- JOBST® Compression Institute serves as a valuable resource for fitters, clinicians, healthcare providers and caregivers interested in learning more about compression therapy and lymphatic and venous disorders.
- Our product experts, called Medical Education Advisors, and our devoted team of Sales Professionals, provide customized education and solutions for both fitters and therapists.



Social Responsibility

- Every year, JOBST® donates hundreds of compression garments to healthcare organizations who share our mission to help improve patient quality of life. Whether it is funding education through organizations such as the Lymphatic Education and Research Network (LE&RN) or helping those recover after a devastating natural disaster, JOBST® does not hesitate to lend a helping hand.

Advanced Wound Care Products

Cutimed offers a wide range of advanced wound care products for a variety of wound phases and exudate levels that result in effective wound healing.

Cutimed® Siltec®

For managing moderate to high levels of exudate

Cutimed® Siltec® has been designed for non-interrupted wound healing and safe fluid handling.

- Gentle adhesive silicone wound contact layer provides protection of fragile new tissue
- Absorbs and retain fluid even under compression



Cutimed® Siltec® Sorbact®

For helping to manage infection and removing exudate

Cutimed® Siltec® Sorbact® is an innovative combination of absorption and infection management. The super-absorbent particles within the foam layer effectively adjust to absorb moderate to high, viscous exudate and retain it even under compression.

- Increased absorption due to super-absorbent hot-melt stripes



Cutimed® Sorbion®

For managing moderate to high levels of exudate

Cutimed® Sorbion® is a line of super-absorbent dressings that absorbs and retains exudate to support a moist healing environment, even under pressure.

- Sorbion Sana helps protect healthy and sensitive tissue
- Sorbion Sachet® provides soft debridement for wound bed preparation



Cutimed® Sorbion® Sorbact®

For removing bacteria and high levels of exudate

Cutimed® Sorbion® Sorbact® is a hydrophobic microbe binding dressing that safely utilizes a physical mode action versus a chemical mode of action.

- The super-absorbent core absorbs and retains moderate to high levels of exudate
- Reduces the risk of maceration even under compression



Compression Bandages

A complete range of bandage products for the management of lymphedema. Designed with the patient and therapist in mind for optimal Complete Decongestive Therapy (CDT) outcomes.

Multi-Layer Bandaging



Comprilan®

- Short stretch compression bandage
- Provides high working pressure and low resting pressure
 - Easy care and durable



CompriFoam®

- Open cell foam bandage
- 100% polyurethane foam
 - Reusable



Tricofix®

- Lightweight stockinette
- For use under Plaster of Paris or synthetic casts
 - Breathable



Artiflex®

- Non-woven padding bandage
- Used in Lymphology bandaging
 - Provides padding around bony or irregular shapes
 - Conformable and air-permeable
 - Supports even distribution of pressure



Elastomull®

- Soft comfortable bandage
- Ideal for fingers and toes
 - Available in sterile or non-sterile



Isoband®

- Elastic multipurpose bandage
- Short stretch with longitudinal elasticity for high working pressure
 - Reusable
 - Porous weave provides good ventilation to the skin

JOBST® Compri2 2-Layer Bandaging

The JOBST® Compri2 is a 2-Layer Bandage System* specifically designed for easy application. It provides graduated and reliable compression for up to seven days to help achieve effective edema reduction and helps promote wound healing.



1st Layer - Non Woven Padding

- **Stretchable** - simplifies application
- **Active Moisture Management** - helps keep skin dry and comfortable
- **Customized Padding** - allows for anatomical adjustments and additional protection to sensitive areas



2nd Layer - Short Stretch Bandage

- **Cohesive** - delivers reliable compression
- **Adjustable** - allows modification during application
- **Visual** - indicator offers simple and effective application support

*As JOBST® Compri2 does not include a wound dressing, any wounds must be treated with a suitable wound dressing before the compression bandage is applied.

JOBST® FarrowWrap®

Adaptable and Affordable, Comfortable and Containing

Product Range

LITE

A double-laminated fabric suitable for mild to moderate compression needs, making it especially suitable for home care use / self-management patients. The fabric is durable.

Available versions

- Upper and lower limb
- Ready-to-Wear (RTW)
- Trim-To-Fit (TTF)
- Custom-Made (CM)

20-30 mmHg* Mild to moderate chronic edema / lymphedema




STRONG

Premium fabric with strong compression and soft inner layer suitable for patients with a good limb shape and moderate to severe edema.

Available versions

- Upper and lower limb
- Ready-to-Wear (RTW)
- Trim-To-Fit (TTF)
- Custom-Made (CM)

30-40 mmHg* Moderate to severe chronic edema / lymphedema with a good limb shape




CLASSIC

The original material that has been proven over many years of use on thousands of patients worldwide.

Available versions

- Upper and lower limb
- Ready-to-Wear (RTW)
- Custom-Made (CM)

30-40 mmHg* Moderate to severe chronic edema / lymphedema



BASIC

A Value fabric with quality compression. Uses double-sided VELCRO® Brand product so it is possible to change band length and trim as limb reduces.

Available versions

- Lower limb
- Ready-to-Wear (RTW)

30-40 mmHg* Moderate to severe chronic edema / lymphedema




4000

The inner sleeve and 4-band design make for easy application. Anatomically contoured to help prevent gapping.

Available versions

- Lower limb
- Ready-to-Wear (RTW)

30-40 mmHg* Moderate to severe chronic edema / lymphedema



JOBST® Elvarex®

JOBST's Gold Standard



JOBST® Elvarex® represents JOBST's gold standard in custom daytime garments for the management of lymphedema. These flat-knit garments offer the very best in comfort and containment for even the most severe cases. With a natural rubber core*, Elvarex® provides unsurpassed structure and control for patients of all shapes and sizes.

JOBST® Elvarex® - Key Features

- Flat-knit structure bridges skin folds instead of working into them
- Natural rubber latex core provides superior containment and allows garment to return to original shape for extended wear
- Excellent breathability for comfort and patient compliance
- Seven compression classes for optimal therapeutic outcomes
- Vast product assortment covers a diverse range of specifications and requirements of patients
- SPF 80 UV protection

















Basic Styles



Compression classes

Class 1	Class 2
Class 2 forte	Class 3
Class 3 forte	Class 4
Class 4 super	

Extremities / Colours

Upper	1-2	 Beige	 Black	 Navy	 Honey	 Cranberry	 Cherry	 Caramel		
	2-Forte	 Beige	 Black	 Navy	 Honey	 Cranberry	 Cherry			
Lower	1-4 Super	 Beige	 Black	 Navy	 Cocoa	 Cranberry	 Grey	 Henna**	 Graphite**	 Denim**

*CAUTION: These products contain natural rubber latex which may cause allergic reactions.

** Not available in Toe Cap

JOBST® Elvarex® Plus*


The Cutting Edge of Flat-Knit Technology

JOBST® Elvarex® Plus uses advanced flat-knitting technology to construct garments in three dimensions to more closely conform to the shape of patients' hands and feet. These seamless garments help minimize potential irritation between the fingers and toes, and at the base of the thumb, help patients with the comfort they want and the performance they need.


*Caution: This product contains natural rubber latex which may cause allergic reactions.



Basic Styles



Compression classes



Type / Colours

Hand	Beige	Black	Navy	Honey	Cranberry	Red	Caramel	Cherry
Foot	Beige	Black	Navy	Cocoa	Cranberry	Red	Grey	

*CAUTION: These products contain natural rubber latex which may cause allergic reactions.


JOBST® Elvarex® Soft

Gentle and Not Made with Natural Rubber Latex


JOBST® Elvarex® Soft is the solution for patients wanting a softer, gentler compression garment. The flat-knit fabric allows for easier donning and meets the specific needs of many elderly and pediatric patients, or patients with lipedema. While soft and forgiving, Elvarex® Soft still provides the containment that mild to moderate lymphedema patients need.



Basic Styles



Compression classes



Extremities / Colours

Upper CCL 1-2	Beige	Black	Navy	Cranberry	Cherry	Honey	
Lower CCL 1	Beige	Black	Navy	Cranberry	Cherry	Grey	Cocoa

JOBST® Confidence



Next Generation of Made-to-measure Flat-knit Compression Garments

JOBST® Confidence helps your patients feel unrestricted, with a one-of-a-kind fit garment that's as unique as they are. The new generation of flat-knit technology used in JOBST® Confidence allows for a precise fit to the real anatomical body shapes. This personally refined fit not only provides comfort, but also supports your patients with a greater degree of movement.

JOBST® Confidence - Key Features

- 3D innovative knitting method allows the garment to conform to individual body shapes
- Conforming fabric that also offers firm support
- Moisture wicking with double-layer construction for optimal comfort
- Flexible, conforming fabric to make putting on and taking off easier

Basic Styles



Compression classes


Upper extremity

Class 1
Class 2


Lower extremity

Class 1
Class 2
Class 2

New Design Options



Deco Line
The deco line is positioned on the front of the garment



Patient Initials
Initials are located under the top band of the garment

Comfort Options



Functional zone*



Slant toe



Hallux valgus closed toe



Hallux valgus open toe

*Available for knee, ankle and elbow

Colours

Beige	Red Heather
Black	Anthracite Heather
Caramel	Jeans Heather

JOBST® Elvarex® Soft Seamless

Soft and Seamless, Where it Matters Most

JOBST® Elvarex® Soft Seamless combines all of the benefits of Elvarex® Soft with the seamless construction that patients want in a glove or toe cap. Not only are these flat-knit garments soft and not made with natural rubber latex, they are seamless, improving comfort by helping to eliminate potential irritation between the fingers and toes.

JOBST® Elvarex® Soft Seamless - Key Features

- Soft, durable, and breathable flat-knitted fabric
- Flat edge finish at fingers and toes decreases bulk and prevents curling
- Easy to don, comfortable to wear

Basic Styles



Compression classes

Class 1	Class 2
---------	---------



Extremities / Colours

Upper CCL 1-2	Beige	Black	Navy	Cocoa	Honey	Cranberry	Cherry
Lower CCL 1-3	Beige	Black	Navy	Cocoa	Grey	Cranberry	Cherry

JOBST® Custom Seamless Soft

Custom-made Compression Hosiery

JOBST® Custom Seamless Soft is a custom-made circular knit garment that looks and feels like hosiery. Appropriate for mild to moderate edema and vascular disorders, these garments wick away moisture and help keep skin cool and dry. With Custom Seamless Soft, patients don't have to compromise on fit or style to wear medical grade compression.


JOBST® Custom Seamless Soft - Key Features

- Multi-fiber double covered inlay yarns increase softness and improve donning
- Refined opacity covers blemishes and scars
- High stretch core yarns increase recoil (stretch and recovery) to prevent sliding
- Reinforced foot zones increase strength and durability with added comfort


Colours

 Sand	 Black	 Espresso	 Caramel
 Navy	 Bronze	 Cranberry	

Basic Styles



Compression classes



JOBST® Bellavar®

Firm Circular Knit Containment

JOBST® Bellavar® is designed to manage moderate to severe edema in the legs by providing controlled stretch and resistance. A custom-made circular knit garment, Bellavar® is seamless and provides a greater level of containment than JOBST® Custom Seamless Soft.


JOBST® Bellavar® - Key Features

- Designed to reduce wrinkling behind the knee
- Variety of styles and options to satisfy patient needs


Colours

 Sand	 Black	 Sun Bronze
--	---	--

Basic Styles



Compression classes



JOBST® Relax

Redefining Lymphedema Night Care

JOBST® Relax is the revolutionary nighttime garment that is designed for more than just sleeping. Its spacer fabric technology is unique in the marketplace and creates a supremely lightweight, flexible, and slim-line garment. Relax brings exceptional quality of life benefits and style to patients with mild to moderate lymphedema through the use of advanced knitting techniques and CoolMax® yarns. Because Relax is so breathable and flexible, patients can remain compliant through the evening hours.



JOBST® Relax - Key Features

- Custom-made, flat-knit technology provides gradient compression therapy and stays in place
- Textured surface and monofilament spacer combine to create differential pressure areas that support stimulation of lymph flow while a patient is at rest
- The integrated zipper option with stay-down puller helps facilitate donning and doffing. The fabric underneath the zipper protects the skin

Arm Styles




Armsleeve
(CG1)



Armsleeve and
Gauntlet Combined
(AG1)

Leg Styles



Knee
(AD)



Thigh
(AG)

Colours



Compression classes

15-20 mmHg*

20-30 mmHg*

JOBST® Custom Seamed

Nearly Limitless Solutions for Unique Challenges

JOBST® Custom Seamed custom-made garments were the first of their kind to manage lymphedema, and while much has changed in the past 70 years, our ability to find solutions to unique challenges has not. Custom Seamed garments are available for every part of the body and in countless configurations. Every garment is precisely cut and sewn by hand in Hickory, North Carolina.

JOBST® Custom Seamed - Key Features

- Durable fabric for long lasting wear
- Precise gradient compression based on computer design
- Variety of options: zippers, linings, silicone bands and more

Colours



Compression classes

15-65+

Basic Styles



Armsleeve



Glove



Gauntlet



Knee



Thigh



Waist High



Capri



Bermuda



Face Mask



Sleeveless
Vest



Body Brief
with Sleeves



Abdominal
Pelvic Garment

JOBST® Bella™ Lite


The Entry Point to JOBST® Quality

JOBST® Bella™ Lite upper extremity garments are designed for patients with mild to moderate lymphedema who are recently diagnosed, post-operative, or do not need the containment of a flat-knit garment. The lightweight fabric is breathable, designed to be comfortable, and easy to don - perfect for patients learning to make garment compliance part of their daily routine.

JOBST® Bella™ Lite - Key Features

- Fits a wide range of arm and hand shapes and available with or without a silicone band
- Two armsleeve lengths for optimal fit
- Fine seams with ultra thin and smooth compression fabric

Basic Styles




Armsleeve

Gauntlet

Colours



* Not available in all compression levels and styles



Compression classes

15-20 mmHg*

20-30 mmHg*

JOBST® Bella™ Strong

Comfortable, Cost-effective Compression

JOBST® Bella™ Strong upper extremity garments are the ready-to-wear solution for patients who need to manage their lymphedema, but do not need the containment provided by a flat-knit custom garment. Bella™ Strong makes JOBST® quality available to a large number of patients by offering a wide range of sizes, lengths, different compression classes, and more accessible pricing.

JOBST® Bella™ Strong - Key Features

- Covers wide range of patient dimensions
- Designed to be easy to don, armsleeve comes with donning slip
- Comfort flat wrist seam reduces rolling
- Cotton-like feel

Basic Styles



Armsleeve

Gauntlet

Colour



Compression classes

15-20 mmHg*

20-30 mmHg*

30-40 mmHg*

* All compressions shown are the mean compression for an average forearm size (armsleeve) or mean compression for an average palm size (glove / gauntlet). Not all styles may be available in all compression ranges.

JOBST® SoftFit™

Our innovative JOBST® SoftFit™ technology is a unique system of silicone yarns knitted into the welt of our compression garments. It provides gentle, seamless friction on the skin, decreasing the chance of constriction and slipping down.



LOOKS GOOD, FEELS GOOD

With the SoftFit™ band, garments are not only visually appealing and more comfortable to wear, they provide precisely the right amount of pressure, right where it is needed.

IMPROVED COMPLIANCE, IMPROVED QUALITY OF LIFE

Compression garments that feel comfortable lead to greater patient compliance.

MORE OPTIONS, HAPPIER PATIENTS

SoftFit Technology is offered in a wide variety of compression garments and styles.

Band Variations

Product	Dotted Silicone	Lace Silicone	SoftFit	Sensitive	Ribbed Elastic Band
Elvarex® Upper Extremity	2.5 cm, on top, inside 5 cm, on top, inside		3 cm, on top (CCL 1 and 2 only)		
Elvarex® Lower Extremity	2.5 cm, on top, inside 5 cm on top, inside 5 cm on top, waistband		5 cm, on top (Knee High only, CCL 1 thru 3)	5 cm, on top (Waist High only)	5cm, on top (Waist High only)
Elvarex® Soft Upper Extremity	2.5 cm, on top, inside 5cm, on top only		3 cm, on top		
Elvarex® Soft Lower Extremity	2.5 cm, on top (Knee High only) 5 cm, on top		5 cm, on top (Knee High only)		
Confidence	5 cm, on top (Thigh High only) 2.5 cm on top for Armsleeve		5 cm, on top (Knee High only) 2.5 cm on top for Armsleeve		
Bella™ Lite Armsleeve	5 cm, on top				
Bella™ Strong Armsleeve	5 cm, on top				
Custom Bellavar	5 cm, on top (Knee & Thigh High) 2.5 cm, on top (Knee High only)	6 cm, on top (Thigh High only)			
Custom Seamless Soft	5 cm, on top (Knee & Thigh High) 2.5 cm, on top (Knee High only)	6 cm, on top (Thigh High only)	5 cm, on top, (Knee High only) (CCL 1 and 2 only)	5 cm, on top, (Thigh High only)	
Custom Seamed	2.5 cm, on top 5 cm, on top			5 cm, on top	

Functional Zones for JOBST® Elvarex® & Elvarex® Soft

Unique and innovative Functional Zones are created by the use of technologically advanced knitting techniques and specialized yarn selection. Functional Zones increase patient comfort through increased flexibility and breathability, without compromising medical effectiveness.

JOBST® Elvarex® Functional Zones

- Use of highly specialized yarn and knitting techniques.
- More comfortable wearing due to increased breathability and flexibility.
- Special zones combat wrinkling when bending and speed moisture dissipation.

Functional Zone "Pressure Panty Elvarex® (body-bandage)"



- Available in:
- Body bandage: CCL 1-3 Forte
 - Corresponding leg stockings available in all CCL (1-4 Super)

Functional Zone "Top Comfort Zone"



- Available in:
- B'G, BG, CG, AG, AT Pressure panty Elvarex® (body bandage), Bermuda, Capri pants, stump stockings (ST)

Functional Zone "Knee Comfort Zone"



- Available in:
- CCL 2-3

Functional Zone "T-Heel"



- Available in:
- CCL 2-3 Forte

Functional Zone "Slant Foot"



Functional Zone "Bias top end"



- Available in:
- All compression classes (CCL 1-2 Forte)

Functional Zone "Elbow Comfort Zone"



- Available in:
- CCL 2
 - Exception: not available for "single piece hand/arm"

JOBST® Donning Aids



JOBST® Donning Gloves*

Putting on JOBST® compression garments is always easier with JOBST® Donning Gloves! They help create friction, provide grip and help to evenly distribute the fabric on the leg.

* Caution: This product contains natural rubber latex which may cause allergic reactions.



JOBST® Stocking Donner

Easy to use device to help apply JOBST® compression garments. Ideal for people who have limited dexterity or have difficulty reaching their feet. Use with closed and open toe styles. Available in 2 sizes: Regular and Large Calf.



Arion Application Aids

Arion Application Aids assist in the donning and doffing of upper and lower extremity garments. These application aids can be used for both open ended and closed ended garments, and can be very helpful for patients with tenderness or weakness in their donning limb or are not strong enough to don the garment themselves.



Stocking Slip

For use with open toe styles only. Helps the stocking slide smoothly onto the foot by reducing friction between the patient's skin and the compression stocking.



JOBST® Roll-On Adhesive

- Designed to keep JOBST® Hosiery in place
- Ideal for thigh high lengths that do not have a silicone band or when a garter belt is not an option
- Helpful Hint: For flat-knit garments, apply in vertical strips



The JOBST® Compression Institute serves as a resource for healthcare providers and clinicians seeking further education about vein diseases and lymphedema. This website provides literature, videos, and links to our national affiliates. Additionally, the JOBST® Compression Institute offers a variety of nationwide workshops and continuing education classes to train you on JOBST® products and proper measuring and fitting techniques. Up to 9 hours of CEU credits available. Register at <http://www.jobstcompressioninstitute.com>.

PREFER TO SHOP ONLINE?

JOBST® eShop

- Place your JOBST® Custom Lymphology orders online
- Intelligent patient management
- Easy measurement documentation
- Time efficient ordering



To get started, visit <https://eshop.jobst-ca.com/>

Questions?

Email us at ca.elvarex@essity.com or call us at 1-877-978-5526.

Sorbact® is a registered trademark of ABIGO Medical AB.
VELCRO® is a registered trademark of Velcro BVBA.
CoolMax® is a registered trademark of INVISTA.
Easy-Slide® is a registered trademark of Arion International.
Please refer to the product label and / or package insert for full instructions on the safe use of these products.



JOBST®,
an Essity brand

 jobstcanada.com essity.com

Essity
1275 North Service Road West, Suite 800
Oakville, ON Canada L6M 3G4
Tel. 1-877-978-5526 Fax 1-877-978-9703
52209E R17 136342E ©Essity Hygiene and Health AB, 2022 C22