

DEFY GRAVITY

A woman with dark hair pulled back, smiling and posing with one hand on her hip. She is wearing a vibrant blue, form-fitting, off-the-shoulder mini dress. Underneath the dress, she wears grey long-sleeved compression sleeves and grey leggings. She is also wearing white sneakers. The background is a solid, muted blue.

SHAPE THEIR COMFORT

JOBST Confidence®

with a personalized fit, thanks to our unique Contour Fit technology, your patients can be as active as they want to be.

**"WEARING COMPRESSION IS
ENABLING ME TO LIVE ENTIRELY
FREELY AND NOT THINK OF MY LEGS
WHEN I'M OUT AND ABOUT."**

Paola - LIPEDEMA



JOBST CONFIDENCE CAN MAKE A DIFFERENCE

You are passionate about finding the right compression garments for your patients, as you don't want them to be held back by their condition. JOBST Confidence is designed with a clear purpose: to allow your patients to feel comfortable throughout the day, in everything they do. That's why it has Contour Fit technology and a specific construction for moisture management. These garments are ideal for patients with mild to moderate lymphedema and lipedema*. They precisely fit the individual anatomical shapes, giving them greater freedom of movement.

THE TECHNOLOGY BEHIND OUR COMFORT


CONTOUR FIT

Unlike traditional flat-knit technology, Contour Fit works by adapting the fit of the garment at four positions on the limb (instead of one) to create a garment that truly contours the unique form of your patient's body. This results in special Contour Points, that show where the garment has been individually tailored to allow your patients to move more freely.

MOISTURE MANAGEMENT

JOBST Confidence garments are designed to keep the skin dry. This is due to a unique construction and yarn composition.

*Without pronounced skin folds and shape distortions



JOBST CONFIDENCE: FOR MORE FEEL GOOD MOMENTS

Innovative knitting
method allows for high
ability to conform to
individual body shapes



Made with a supple
conforming fabric that also
gives the reassurance of
firm support



Advanced moisture
management system with
double-layer construction for
optimal wearing comfort



UNIQUE BODIES DESERVE A UNIQUE FIT

For the first time ever, JOBST Confidence is able to provide a dedicated knitting pattern for both conditions.

WHAT'S YOUR PATIENTS' SHAPE?

LIPHEDEMA



LYMPHEDEMA



Measuring

Measure **both** body types according to our guideline.



Recognizing

The manufacturing system interprets measurements and recognizes the body shape.



Producing

After the knitting pattern selection the innovative **Contour Fit** knitting technique produces a garment with an even **better** personalized fit.

A CLINICAL STUDY SHOWS THE SUPERIOR COMFORT OF JOBST CONFIDENCE*

In a clinical study the features below were rated by patients in terms of their satisfaction at the beginning of the study with regards to their previously worn product and at the end of the study with regards to JOBST Confidence. These features were evaluated better for JOBST Confidence and the number shown display the rating for JOBST Confidence.

83.5%

COMFORT

very good or good*

80%

FIT

very good or good*

91.8%

MOISTURE MANAGEMENT

very good or good*

84.7%

FREEDOM OF MOVEMENT

very good or good*

*Prospective multicentre study to evaluate the clinical performance and safety of innovative, made-to-measure compression garments in daily routine. n=85. In preparation for publication, 2023

JOBST CONFIDENCE - DESIGNED FOR SUPERIOR COMFORT

Having a comfortable garment is vital as it means that your patient is more likely to wear their garment regularly, which is important for effective treatment. JOBST Confidence offers besides its unique design, options that are also designed to make it more comfortable.



PANTY DESIGN FOR COMFORTABLE FIT

Our panty styles feature an ergonomic seam placement that better follows the natural shape of lymphedema and lipedema patients' bodies. This can help reduce discomfort and restriction in their daily lives.

A DEDICATED DESIGN FOR MEN

All panty styles are also available with an ergonomic seam placement that follows the natural shape of men's bodies. The fly option, specifically designed for men, has an improved fit and is designed for higher wearing comfort.



COMFORT OPTIONS

Functional Zones are designed to enable greater mobility around the elbow, knee and ankle.



Functional zone knee*



Functional zone arm*



Slant toe



Hallux valgus
closed toe




Hallux valgus
open toe


*Available for knee, ankle and elbow

JOBST CONFIDENCE - DESIGNED FOR INDIVIDUAL NEEDS

DESIGN OPTIONS



Deco line



Initials



Floral
waistband (5 cm)

COLOURS




Beige



Caramel



Black



Jeans Heather



Anthracite Heather



Red Heather

STYLES



AB1 Sock



AD Knee high



AG Thigh high



CG Arm sleeve



CG1 arm sleeve
with gauntlet



AT Panty



BT Legging



BIT Legging



CT Capri



AG-HT One
legged tights



ET Bermuda

COMPRESSION CLASSES

Upper limb

CCL 1 ▶ 15 - 21 mmHg

CCL 2 ▶ 23 - 32 mmHg

Lower limb

CCL 1 ▶ 18 - 21 mmHg

CCL 2 ▶ 23 - 32 mmHg

CCL 3 ▶ 34 - 46 mmHg



JOBST®,
an Essity brand



Essity
1275 North Service Road West, Suite 800
Oakville, ON Canada L6M 3G4
Tel. 1-877-978-5526 Fax 1-877-978-9703

64037E R3 136406 ©Essity Hygiene and Health AB, 2023 K23