



FarrowWrap® Upper Extremity

Ready-To-Wear • Custom • Trim-to-Fit

TO ORDER:

Tel: 1-877-978-5526

1-877-358-2739

Fax: 1-877-978-9703

Patient Name _____ DOB _____ Date _____

Address _____ Gender M ☐ F ☐

City/Province/Postal Code _____

Diagnosis _____

Doctor / Address _____

City/Province/Postal Code _____

PO#

Original Order ☐ Reorder w Changes ☐
Exact Reorder ☐

Ready-To-Wear Arm

Fabric

☐ LITE

Colour

☐ Beige (00)

☐ Black (01)

Length

☐ Short (S)

☐ Regular (R)

☐ Long (L)

Size

☐ Small (2)

☐ Medium (3)

☐ Large (4)

Side

☐ Left (L)

☐ Right (R)

Custom

Style

☐ Custom

Fabric

☐ Classic

☐ LITE

☐ STRONG

Colour

☐ Tan

Side

☐ Left (L)

☐ Right (R)

Trim-To-Fit

Size


☐ Medium (3)

☐ Extra Large (5)

Side

☐ Left (L)

☐ Right (R)



Circumference (c)	Left	Right
cG (axilla)*		
cF ¹ (widest bicep)		
cF (distal bicep)		
cE (elbow crease)*		
cD ¹ (widest forearm)		
cD (distal forearm)		
cC (wrist)*		

Lengths	Left	Right
IC-G (wrist to axilla)*		
IC-E (wrist to elbow)		

Comments

*Circumference measurements needed for Ready to Wear garments.
Please note that all measurements are needed for Custom garments.



FarrowWrap® Upper Extremity

Ready-To-Wear • Custom • Trim-to-Fit

READY-TO-WEAR ARM

Size	Wrist (cC)	Elbow (cE)	Axilla (cG)
S	14-18 cm	20-27 cm	22-31 cm
M	16-21 cm	25-34 cm	29-39 cm
L	19-25 cm	30-40 cm	32-45 cm
Length (C-G)	Short 40-43 cm		
	Regular 44-47 cm		
	Long 48-51 cm		

TRIM-TO-FIT ARM

Size	Wrist (cC)	Forearm (cD)	Elbow (cE)	Bicep (cF)	Axilla (cG)	Length (C-E)	Length (E-G)
M	15-31 cm	20-36 cm	20-38 cm	25-42 cm	25-44 cm	20-25 cm	15-23 cm
XL	15-35 cm	20-40 cm	25-42 cm	25-45 cm	30-48 cm	20-28 cm	

READY-TO-WEAR ARM - TAN

Code	Description	Quantity
7665035	LITE Armpiece Long Tan Small Left	
7665000	LITE Armpiece Long Tan Small Right	
7665001	LITE Armpiece Long Tan Medium Left	
7665002	LITE Armpiece Long Tan Medium Right	
7665003	LITE Armpiece Long Tan Large Left	
7665004	LITE Armpiece Long Tan Large Right	
7665005	LITE Armpiece Regular Tan Small Left	
7665006	LITE Armpiece Regular Tan Small Right	
7665007	LITE Armpiece Regular Tan Medium Left	
7665008	LITE Armpiece Regular Tan Medium Right	
7665009	LITE Armpiece Regular Tan Large Left	
7665010	LITE Armpiece Regular Tan Large Right	
7665011	LITE Armpiece Short Tan Small Left	
7665012	LITE Armpiece Short Tan Small Right	
7665013	LITE Armpiece Short Tan Medium Left	
7665014	LITE Armpiece Short Tan Medium Right	
7665015	LITE Armpiece Short Tan Large Left	
7665016	LITE Armpiece Short Tan Large Right	

READY-TO-WEAR ARM - BLACK

Code	Description	Quantity
7665023	LITE Armpiece Regular Black Small Left	
7665024	LITE Armpiece Regular Black Small Right	
7665025	LITE Armpiece Regular Black Medium Left	
7665026	LITE Armpiece Regular Black Medium Right	
7665027	LITE Armpiece Regular Black Large Left	
7665028	LITE Armpiece Regular Black Large Right	
7665029	LITE Armpiece Short Black Small Left	
7665030	LITE Armpiece Short Black Small Right	
7665031	LITE Armpiece Short Black Medium Left	
7665032	LITE Armpiece Short Black Medium Right	
7665033	LITE Armpiece Short Black Large Left	
7665034	LITE Armpiece Short Black Large Right	
7665017	LITE Armpiece Long Black Small Left	
7665018	LITE Armpiece Long Black Small Right	
7665019	LITE Armpiece Long Black Medium Left	
7665020	LITE Armpiece Long Black Medium Right	
7665021	LITE Armpiece Long Black Large Left	
7665022	LITE Armpiece Long Black Large Right	

CUSTOM CLASSIC

Code	Description	Quantity
7599500	Classic Arm	

CUSTOM STRONG

Code	Description	Quantity
7598700	STRONG Arm	

CUSTOM LITE

Code	Description	Quantity
7599100	LITE Arm	

FARROWWRAP TRIM-TO-FIT

Code	Description	Quantity
7612000	LITE Armpiece (Left) Medium with 2 Standard arm liners	
7612002	LITE Armpiece (Right) Medium with 2 Standard arm liners	
7612001	LITE Armpiece (Left) Extra Large with 2 Standard arm liners	
7612003	LITE Armpiece (Right) Extra Large with 2 Standard arm liners	



JOBST®,
an Essity brand



Essity
1275 North Service Road West, Suite 800,
Oakville, ON L6M 3G4
Tel 1-877-978-5526 Fax 1-877-978-9703
62363E R1 136122E ©Essity Hygiene and Health AB 2022 H22